

FAHL® PERCUTAN DILATION SET

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de

PICTOGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

Catalogue number

Batch code

Use by

Date of manufacture

Consult instructions for use

Keep away from sunlight

Keep dry

Do not reuse

Medical devices

Sterilized using ethylene oxid

Do not resterilize

Do not use if package is damaged

Simple sterile barrier system

CE marking identification number of the notified body

Manufacturer

FAHL® PERCUTAN DILATION SET

GENERAL INFORMATION

These instructions are applicable for the **FAHL®** Percutan Dilation Set.

The user manual serves to inform users in order to ensure technically correct handling.

Read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Keep the package for as long as you use the product.

It contains important information on the product!

1. INTENDED PURPOSE

FAHL® Percutan Dilation Sets are kits for placement of a tracheostoma by percutaneous dilation tracheostomy.

2. SAFETY INFORMATION

Percutaneous dilation tracheostomy is a procedure that must be performed by doctors only or under the instruction and supervision of doctors who are trained in and have experience with using the Ciaglia method and/or the Griggs method.

At least two surgeons must carry out the procedure:

1. To perform the operation

2. For monitoring with the bronchoscope and to monitor the anaesthesia and vital signs of the patient.

The product is intended for use in a clinical environment.

Before starting the procedure, identify the tracheal and pretracheal structures with ultrasound or similar to determine the optimum puncture site.

Monitor the intratracheal position of the puncture cannula, the guide wire, the dilator or the tracheal spreader forceps and the tracheal cannula to minimise the danger of paratracheal insertion or injury to the trachea.

Use an endotracheal tube during the procedure to maintain the patient's respiration with continuous oximetric monitoring.

Take precautions in case of complications to ensure that the ventilation for the patient can be maintained through alternative airways (e.g. intubation).

Have a spare **FAHL®** Percutan Dilation Set and a spare tracheal cannula ready in case of complications.

If haemorrhage occurs, aspirate the blood immediately to prevent it from accumulating in the airways.

The product is a single-patient product and is intended solely for single use.

Follow the instructions for use when used in combination with other medical devices.

3. INDICATION

- Long-term ventilation
- Airway obstruction

- Injuries in the region of the larynx or mouth
- Tumours in the region of the larynx or mouth

4. CONTRAINDICATION

Absolute contraindications:

- Respiratory emergencies
- Age < 18 years
- Unstable cervical fractures
- Infections in the tracheostomy region
- Malignancies in the tracheostomy region
- Unclear anatomical conditions
- Difficult intubation conditions

Relative contraindications:

- Coagulation disorders
- Oxygenation disorders
- Enlarged thyroid
- Prior surgical procedures in the neck region (e.g. thyroidectomy)
- Haemorrhagic diathesis (increased tendency to haemorrhage)
- Deep-seated airways
- Immunosuppression

5. COMPLICATIONS AND ADVERSE REACTIONS

- Venous haemorrhage
- Tracheal injuries
- Rupture of the tracheal cartilage
- Oxygen saturation
- Pleural effusion
- Infections of the stoma
- Herniation
- Tracheo-oesophageal fistula

6. PRODUCT DESCRIPTION

REF 19980: The dilation set consists of scalpel (**A**), syringe (**B**), puncture cannula (14 gauge) with catheter (**C**), guide wire (**D**), pre-dilator (**E**), guide catheter (**F**), dilator with hydrophilic coating (**H**), 4 compresses (**G**) and sterile wipe (**I**).

REF 19981: The dilation set consists of scalpel (**A**), syringe (**B**), puncture cannula (14 gauge) with catheter (**C**), guide wire (**D**), pre-dilator (**E**), dilation forceps (**J**), 4 compresses (**G**) and sterile wipe (**I**).

REF 19982: The dilation set consists of scalpel (**A**), syringe (**B**), puncture cannula (14 gauge) with catheter (**C**), guide wire (**D**), pre-dilator (**E**), guide catheter (**F**), dilator with hydrophilic coating (**H**), obturator sizes 7, 8 and 9 (**K**) and 4 compresses (**G**) and sterile wipe (**I**).

REF 19920, 19921, 19922, 19923: FAHL® Percutan Dilation Set (**19980**) with Spiraflex P cannula in various sizes

REF 19990: The guide wire (**D**) has a lock and an insertion guide.

REF 19991: The guide wire (**D**) has a lock and an insertion guide.

The guide catheter (**F**) has a safety lock.

7. APPLICATION

7.1 PREPARATION

The product is intended for use in an aseptic environment.

Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged.

Do not use the product if the packaging has been damaged.

Check the use-by or expiry date.

Do not use the product after this date.

Before use check that the product is complete, intact and in functional condition.

Damaged products must not be used and must be discarded immediately.

The surgeon should select the tracheal cannula appropriate for the procedure.

Prepare the tracheal cannula as directed by the instructions for use.

Open the packaging.

Check that the guide wire can be easily pushed through the catheter of the puncture cannula, the guide catheter and the pre-dilator.

REF 19980, 19982: Check whether the tip of the dilator with hydrophilic coating is blocked by the safety lock of the guide catheter.

CAUTION!

The tip of the dilator with hydrophilic coating must not be pushed past the safety lock of the guide catheter!

REF 19981: Check that the guide wire can be guided through the tracheal spreader forceps without difficulty.

Dispense sterile water or saline solution into the trough of the blister tray.

Partially fill the syringe with sterile water or sterile saline solution and place it on the puncture cannula with Teflon catheter to prepare the aspiration test.

REF 19980, 19982: Immerse the narrow end of the dilator with hydrophilic coating in the sterile water or saline solution from the tip to the marking to activate the antifriction coating.

OR

Use the included compresses to moisten the front section of the dilator with hydrophilic coating.

CAUTION!

Always keep a bronchoscope at hand for regular checks.

7.2 PREPARATION OF THE PATIENT

1. The patient must be optimally pre-oxygenated immediately before the procedure.
2. Examine the patient for anatomical anomalies that could have a negative effect on the procedure.
3. Position the patient in the supine position with a pillow between the shoulders and head and neck extended (see Fig. 1).
4. Raise the head section of the patient bed **30-40°**.
5. Disinfect the skin in the region of the neck and cover it with sterile cloths.
6. Sedate the patient and monitor in accordance with the hospital guidelines.
7. Aspirate the pharynx before starting the procedure.
8. If an endotracheal tube is used for ventilation, unblock it and pull it back to the entrance to the larynx with the cuff (approx. 1 cm below).

This will prevent perforation of the cuff during the tracheal puncture.

It may be necessary to block the cuff again and adjust the breath volume during ventilation.

7.3 PROCEDURE

WARNING!

The patient must be ventilated during the procedure with an endotracheal tube under continuous oxymetric monitoring.

The puncture of the airways is controlled by aspiration of air.

Tilt the tip of the puncture cannula at an angle of approx. **60°** to the anterior wall of the trachea after the puncture.

This will reduce the risk of injury to the posterior tracheal wall.

Insert the subsequent guide wire and dilator in the same way.

CAUTION!

Make sure that the distal tip of the guide catheter never extends beyond the tip of the guide wire.

The proximal mark of the guide wire must always be at the proximal end of the guide catheter.

This will prevent injury to the posterior tracheal wall.

The safety lock of the guide catheter must be positioned immediately before the tip of the dilator.

Use the marks on the guide catheter to check that the insertion depth is correct.

CAUTION!

Do not insert the dilator past the "MAXIMUM INSERTION" mark (skin level). This will minimise the risk of injury to the trachea at the height of the carina of the trachea.

When inserting the dilator or the tracheal cannula, avoid excessive rotational motions to avoid injuring the trachea.

Perform the dilation in accordance with the selected cannula size so the tracheal cannula can be inserted without requiring significant force.

Obturators:

The insertion aid serves to stabilise the tracheostomy tube during insertion into the tracheostoma.

If you are using a third-party tube, you can use the included obturators (K) to insert it.

Select the obturator to fit the size of the tube.

You should first push the obturator into the tube and then push both over the guide wire into the tracheostoma.

Once the insertion is successfully completed, the obturator can be removed with the guide wire. Make sure that the tracheostomy tube is not pulled out at the same time.

7.3.1 CIAGLA METHOD (19980, 19982)

1. Administer a local anaesthetic to the surgical site.
2. The puncture is generally performed in a caudal direction between the 2nd and 3rd or the 3rd and 4th tracheal ring.
3. Make a 1.5-2 cm incision at the planned insertion point with the scalpel (see Fig. 2) to enable dilation of the skin.
4. Push the puncture cannula with catheter and syringe positioned in a posterior and caudal direction until air bubbles can be aspirated.
5. Tilt the puncture cannula at an angle of approx. 60° to the anterior wall of the trachea (see Fig. 3).
6. Check the position of the puncture cannula in the trachea by aspiration of air in the syringe.
7. Pull the puncture cannula out of the trachea.
Leave the catheter in the puncture cannula.
8. Push the catheter of the puncture cannula several millimetres posterior and caudal (see Fig. 4).
9. Pull the guide wire out of its sleeve (approx. 8 cm) and push it forward until the J-shaped tip is in the straight insertion aid.
10. Position the insertion aid in the catheter of the puncture cannula (see Fig. 5).
11. Attach the dispenser to the catheter of the puncture cannula and push the guide wire through the catheter into the trachea until the second mark is visible in the dispenser window (approx. 10 cm) (see Fig. 5).
12. Make sure that the guide wire can be pushed forward easily and can be freely moved inside the catheter of the puncture cannula.
If the guide wire bends, it must be replaced with a new one.
13. Pull the catheter of the puncture cannula completely out over the guide wire.
Leave the guide wire in position in the trachea.
14. Push the pre-dilator over the guide wire into the trachea.
15. Dilate the puncture site (see Fig. 6).
Make sure that the pre-dilator is aligned and that it does not injure the posterior wall of the trachea.
16. Remove the pre-dilator.
The guide wire remains in position as before.
17. Push the guide catheter with the safety lock forward (in the direction of the arrow) over the guide wire until the safety lock of the guide catheter reaches skin level (see Fig. 7).

CAUTION!

Do not push the safety lock of the guide catheter past the skin level!

Do not insert the dilator any deeper than the "MAXIMUM INSERTION" mark at skin level in order to avoid injury to the trachea and/or the carina.

Avoid rotational motions when inserting the dilator.

18. Push the dilator over the guide catheter until its distal tip reaches the safety lock.

CAUTION!

Do not push the dilator past the safety lock of the guide catheter!

19. Carefully push the dilator along with the guide catheter into the trachea forward and back to stretch the puncture channel.
20. Expand the puncture channel until the opening is slightly larger than the outside diameter of the selected tracheal cannula (see Fig. 8).
Orient yourself with the aid of the marks at 38 Ch/FR und 41 Ch/FR ("MAXIMUM INSERTION") on the dilator. During the dilation ensure that the position of the guide catheter and the guide wire in relation to the dilator is maintained.
21. Remove the dilator (see Fig. 9).
Leave the guide catheter and guide wire in the trachea.
Follow the instructions for use of the tracheal cannula during insertion.

7.3.2 GRIGGS METHOD (19981)

1. Administer a local anaesthetic to the surgical site.
2. Make a 1.5-2 cm incision at the insertion point (see Fig. 2) to enable dilation of the skin.
3. Carefully prepare the anterior wall of the trachea.
If necessary, also free the centre line of the trachea from tissue.
4. Palpate the cricoid cartilage.

5. If necessary, inject additional local anaesthetic into the prepared region.
If the isthmus of the thyroid is in the preparation area, carefully push it in the cranial direction.
6. Use the puncture cannula with catheter and syringe attached partially filled with sterile saline solution or sterile water.
7. Insert the puncture cannula (typically between the 2nd and 3rd or 3rd and 4th tracheal ring).
8. Push the puncture cannula to the posterior and caudal direction until air bubbles can be aspirated.
The puncture cannula should now be tilted at an angle of approx. 60° to the anterior wall of the trachea (see Fig. 3).

CAUTION!

Always use the bronchoscope to monitor the procedure to prevent injury to the posterior tracheal wall!

PLEASE NOTE!

Do not puncture the cuff of the endotracheal tube!

9. If air can be easily aspirated, pull the puncture cannula out of the catheter and then push the catheter of the puncture cannula a few millimetres posterior and caudal (see Fig. 4).
10. Aspirate air from the trachea again with the syringe to check the position of the catheter of the puncture cannula.
11. Adjust the guide wire with the dispenser so the J-shaped tip is in the straight section of the dispenser and does not extend beyond it.
12. Attach the dispenser to the catheter of the puncture cannula and push the guide wire through the catheter into the trachea until the second mark is visible in the dispenser window (approx. 10 cm) (see Fig. 5).

CAUTION!

Never use a damaged guide wire!

The guide wire must be easy to move and it must be easy to move in the catheter of the puncture cannula.
If the guide wire is bent, replace it with a new one.

13. Push the guide wire in the caudal direction while monitoring with the bronchoscope.
14. After positioning the guide wire, pull the catheter of the puncture cannula completely out over the guide wire, leaving the guide wire in position.
15. Push the pre-dilator over the guide wire into the trachea.
16. Dilate the puncture site (see Fig. 6).
Make sure that the pre-dilator is not tilted against the guide wire, otherwise the tip may be damaged.
17. Guide the pre-dilator monitored by the bronchoscope so it does not injure the posterior wall of the trachea.
18. Remove the pre-dilator without changing the position of the guide wire.
19. Allow the tracheal spreader forceps to be locked in the closed position so the jaws are closed completely from the tip to the opening.
20. Insert the guide wire through the opening at the front end of the tracheal spreader forceps and allow it to extend through the opening behind (see Fig. 10).
Hold the free end of the guide wire tight.
21. Insert the tracheal spreader forceps at the same angle as before during insertion of the puncture cannula over the guide wire until the resistance of the anterior tracheal wall can be felt (11).
If the tracheal spreader forceps enter as soon as the initial insertion into the trachea, continue with step 28.
22. Open the tracheal spreader forceps gradually with both hands so the tissue is widely spread.
23. Pull the tracheal spreader forceps back in the open position.

CAUTION!

The guide wire must be easily moved in the trachea during this process.

24. Guide the closed tracheal spreader forceps through the tracheal wall.
The resistance must perceptibly reduce when the tracheal wall is punctured.

CAUTION!

Before continuing, always make sure that the tracheal wall has been penetrated!

25. If you do not feel a reduction in resistance while puncturing the tracheal wall, it is possible that only pretracheal structures are being penetrated!
If necessary, repeat step 24!
26. Push the tracheal spreader forceps further forward in the posterior/caudal direction to widen the tracheal wall.
27. Keep the handles of the tracheal spreader forceps in the median position and lift it to the vertical position so the tip of the tracheal spreader forceps penetrates the tracheal wall further and are in the trachea in the longitudinal direction.

28. Open the tracheal spreader forceps with both hands so the tracheal wall is spread wide enough to insert the tracheal cannula (12).
29. Pull the tracheal spreader forceps back in the open position.
30. Follow the instructions for use of the tracheal cannula during insertion.

7.4 POST-OPERATIVE CARE

Cover the tracheostoma with a sterile wound compress as required.

Do not change the tracheal cannula for the first seven days after the percutaneous dilation tracheostomy, because the tracheostoma is very unstable in the first few days.

If it is really necessary to change the tracheal cannula, use only the Seldinger method along with intubation instruments for the change procedure.

8. HYGIENE INSTRUCTIONS

Cleaning, disinfection or (re-)sterilisation as well as reuse can impair the safety and function of the product and are therefore not permissible.

9. STORAGE

Store in a dry environment away from sunlight and/or heat.

10. SERVICE LIFE

The product is intended for single use on one patient.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, **FAHL** GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies both to damage caused to the products themselves as a result thereof and to all consequential damage caused thereby.

When using the product beyond the period of use specified under Section SERVICE LIFE, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, **FAHL** GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of **FAHL** GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all **FAHL** GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTCs) which can be obtained directly from **FAHL** GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

FAHL® is a trademark and brand of **FAHL** GmbH, Cologne registered in Germany and the EU member states.

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