LARYVOX® PLUG **LARYVOX® PLUG 4MM**

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PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

en

	Catalogue number
	Batch code
	Use by
	Content (in pieces)
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	For reuse on a single patient
	Medical devices
	CE marking
	Manufacturer

GENERAL INFORMATION

These instructions for use are valid for **LARYVOX® PLUG & LARYVOX® PLUG 4MM**.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

1. INTENDED PURPOSE

The **LARYVOX® PLUG & PLUG 4MM** sealing plug serves for the transient use with a defective shunt valve to prevent transprosthetic leakage (centrally through the shunt valve).

2. SAFETY INFORMATION

Patients must be trained by medical professionals on how to use the **LARYVOX® PLUG** safely.

The **LARYVOX® PLUG** is not intended for permanent use and is only used for the transient sealing of a defective shunt valve. It is not a permanent replacement for a functioning shunt valve.

Please consult a specialist as soon as possible to have the shunt valve professionally examined and replaced if necessary.

LARYVOX® PLUG is a single-patient product and thus intended exclusively for use in a single patient.

Re-use of any kind in other patients and thus also reconditioning for re-use in other patients is not allowed.

Only clean the **LARYVOX® PLUG** with cleaning agents approved by the manufacturer.

The **LARYVOX® PLUG** must be cleaned after each use according to the cleaning instructions (see Section CLEANING).

The safety strap on the **LARYVOX® PLUG** must not be changed or removed under any circumstances.

Never force the **LARYVOX® PLUG** or insert it too deeply into the shunt valve. There is a risk that the shunt valve may be pushed/displaced into the oesophagus. In this case, a physician must be consulted immediately.

No tracheostoma accessories (e.g. HME filter cassettes, base plates, tracheostomy tubes) may be used while using the **LARYVOX® PLUG**. Use could result in an accidental change of position of the shunt valve/**LARYVOX® PLUG** and cause aspiration (inhalation) into the lungs.

Do not perform any type of adjustment or repair on the **LARYVOX® PLUG**. In case of damage, the product may not be used and must be disposed of immediately.

Please also refer to the instructions for use of the shunt valve you are currently using.

3. INDICATION

The **LARYVOX® PLUG** is used to transient seal a shunt valve to prevent fluid leakage through the shunt valve until the shunt valve has been properly inspected and replaced if necessary.

4. CONTRAINDICATION

Do not use if the patient is allergic to the material.

5. COMPLICATIONS AND ADVERSE REACTIONS

If there are signs of shunt valve leakage or **LARYVOX® PLUG** malfunction, a physician should be consulted as soon as possible to check the sealing plug and the shunt valve.

The **LARYVOX® PLUG** must not be used with a 4mm size shunt valve.

Only use the **LARYVOX® PLUG 4MM** with a 4mm shunt valve.

6. PRODUCT DESCRIPTION

The **LARYVOX® PLUG** consists of a safety strap (2)a and a sealing plug (2)b.

CAUTION!

Using the **LARYVOX® PLUG** in the shunt valve makes speech temporarily impossible.

7. APPLICATION

Check the use-by or expiry date.

Do not use after this date.

Do not use the product if the packaging has been damaged.

Never use the **LARYVOX® PLUG** if it displays visible deficiencies or faults. Contact the manufacturer in such cases.

The **LARYVOX® PLUG** must be thoroughly cleaned before each use, also before first use.

Always clean your hands with soap and water before using the **LARYVOX® PLUG**.

7.1 INSERTION OF THE LARYVOX® PLUG

Start use of the **LARYVOX® PLUG** in front of a mirror under optimal lighting conditions so that the shunt valve is clearly visible. Then insert the **LARYVOX® PLUG** as follows:

Guide the tip of the insertion (1b) rod into the rear opening of the **LARYVOX® PLUG** (3).

Fasten the safety strap to the insertion rod ((4) / (5) a and b).

Hold the insertion rod with the safety strap facing downwards (6) and carefully insert the **LARYVOX® PLUG** into the shunt valve until it is fixated correctly (7).

CAUTION!

Never force the **LARYVOX® PLUG** or insert it too deeply into the shunt valve. There is a risk that the shunt valve may be pushed into the oesophagus.

Check the secure placement of the **LARYVOX® PLUG** by slightly pulling the safety strap (8).

Now loosen the safety strap from the insertion rod (9).

Slowly and carefully remove the insertion rod from the opening of the **LARYVOX® PLUG** (10) while rotating it slightly.

CAUTION!

The insertion rod must be removed with special care to avoid dislocation (displacement) / aspiration (inhalation) of the shunt valve into the lungs.

Medical assistance must be sought immediately in the event of aspiration.

If necessary, secure the safety strap with medical adhesive tape (11).

7.2 Removing the LARYVOX® PLUG

Always clean your hands with soap and water before removing the **LARYVOX® PLUG**.

Remove the **LARYVOX® PLUG** in front of a mirror under optimal lighting conditions so that the shunt valve and the sealing plug are clearly visible. Then remove the **LARYVOX® PLUG** as follows:

Carefully pull the safety strap to remove the **LARYVOX® PLUG** from the shunt valve (12).

8. HYGIENE INSTRUCTIONS

Clean the **LARYVOX® PLUG** after each use.

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.

CAUTION!

Cleaning agent and disinfectant residues on the product can lead to irritations of the mucous membranes or other health impairments.

8.1 Cleaning

WARNING!

Under no circumstances clean the LARYVOX® PLUG with cleaning agents not approved by the manufacturer. Never use aggressive household cleaning agents, highly concentrated alcohol or denture cleaning agents.

This may be very hazardous to health! The LARYVOX® PLUG may also be destroyed or damaged.

Place **LARYVOX® PLUG** to soak in clear, lukewarm water for approx. 2-3 minutes for cleaning, or use a mild, neutral-pH washing lotion.

Then rinse with plenty of clear water and allow to dry in air.

CAUTION!

Do not use any sharp-edged cleaning utensils!

Always discuss the exact procedure with your doctor/physician.

Please note that the personal cleaning schedule, which may also include additional disinfection if necessary, must always be coordinated with the attending physician.

The **LARYVOX® PLUG** should be immersed in boiling water (100 °C/212 °F) for at least 5 minutes but no more than 10 minutes.

Dry the product thoroughly with a clean and lint-free cloth after wet cleaning.

8.2 Chemical Disinfection Instructions

Cold-disinfection of the **LARYVOX® PLUG** with special chemical disinfectants is possible.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

CAUTION!

Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).

Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This could badly damage or even destroy the **LARYVOX® PLUG**.

8.3 Disinfection steps

For this purpose, OPTICIT® tube disinfectant (REF **331170**) should be used in accordance with the manufacturer's instructions.

After disinfection, the **LARYVOX® PLUG** must be rinsed very thoroughly with sterile saline solution (NaCl **0.9%**) and then dried.

Dry the product thoroughly with a clean and lint-free cloth after wet cleaning.

9. STORAGE

Cleaned products that are not currently in use should be stored in a dry environment in a clean plastic bag and protected from dust, sunlight and/or heat.

10. SERVICE LIFE

The **LARYVOX® PLUG** is a single-patient product.

When used and looked after correctly, in particular observing the instructions for use, cleaning, disinfection and maintenance provided in the instructions for use, the **LARYVOX® PLUG** has an expected service life of up to **12 months**.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section **10**, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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Detta kan erhållas direkt från Andreas Fahl Medizintechnik-Vertrieb GmbH.
LARYVOX® är ett i Tyskland och EU inregistrerat varumärke som ägs av Andreas Fahl Medizintechnik-Vertrieb GmbH, Köln.