
FAHL® HAND FIX

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



32470 / 32472 32482 / 32485

PICTOGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

	Catalogue number		Keep dry
	Batch code		Do not reuse
	Use by		Medical devices
	Date of manufacture		Do not use if package is damaged
	Content (in pieces)		CE marking
	Consult instructions for use		Manufacturer
	Keep away from sunlight		

FAHL® HAND FIX

1. FOREWORD

These instructions are valid for the **FAHL® HAND FIX**.
The instructions for use are for the information of doctors and nursing staff in order to ensure technically correct handling.
Please read the instructions for use carefully before using the product for the first time!
Keep them in an easily accessible place for future reference.
Please keep the package for as long as you use the product.
It contains important information on the product!

2. INTENDED PURPOSE

FAHL® HAND FIX is intended for use with patients in acute danger of self-harm after exhausting all other appropriate measures.

3. SAFETY INFORMATION

Always observe the fixing specifications of your workplace and the relevant legislation of your country.
Legal regulations and medical requirements must be strictly observed.
The supervising doctor must ensure that the patient's state of health is sufficient for the application of **FAHL® HAND FIX**.
In principle, fixing the patient must be cleared legally and if at all possible dispensed with.
FAHL® HAND FIX must be used only with beds, fixed chairs or supports that allow attachment in accordance with the instructions for use.
The products must be attached so they cannot slip.
The correct seating of the fix must be regularly checked.
The fix belt and adjustment mechanisms of beds, mattresses, variable-pressure mattresses etc., supports or chairs may influence one another, resulting in damage, loss of function and risk of injury.
The selected attachment points must not have any sharp edges.
The bed rails must be raised to prevent accidents such as strangulation.
Patients in fixes require regular monitoring, which means that professional observation by qualified personnel trained in the application is absolutely essential.
Aggressive, angry and restless patients or patients in danger of aspiration require continuous monitoring.
The fix must be snug but must not affect respiration and blood circulation.
If **FAHL® HAND FIX** is too loose, the safety of the patients may be endangered.

FAHL® HAND FIX may stretch when soaked and result in the fix coming loose.
The product is a single-use item.
Reuse may result in infection and/or loss of function.

4. INDICATION

- Anxiety
- Delirium
- During recovery phases
- Disorientation/dementia

5. CONTRAINDICATION

- Suicide danger
- Acute aggression to others and self
- Highly agitated patients
- Dislocation or fracture of limbs
- Wound sites
- Psychological or physical problems (that may become acute)

6. COMPLICATIONS AND ADVERSE REACTIONS

- Danger of decubitus
- Contractures
- Weight reduction (no spontaneous food intake)
- Obstipation prophylaxis and pneumonia prophylaxis

7. PRODUCT DESCRIPTION

The hand fixes are available in different versions.
All versions have a Velcro fastening and an integrated fastening strap.
The integrated strap allows direct fixing.
REF 32470, 32472 FAHL® HAND FIX JUNIOR
These hand fixes are intended for newborns, babies and small children.
REF 32482, 32485 FAHL® HAND FIX
These hand fixes are intended for youths and adults.

8. APPLICATION

Check the product for damage before use.
Position collars (see Fig. 1) and close with Velcro (see Fig. 2).
REF 32470, 32472 FAHL® HAND FIX JUNIOR
Thread the straps around the joint and through the eyelet so the Velcro is covered (see Fig. 3).
Then fasten the straps to a suitable support (see Fig. 4).

9. HYGIENE INSTRUCTIONS

Cleaning, disinfection or (re-)sterilisation as well as reuse can impair the safety and function of the product and are therefore not permissible!

10. STORAGE

Store in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

Change when dirty, defective and according to hygiene guidelines.

12. DISPOSAL





The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas **FAHL** Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.
In particular, Andreas **FAHL** Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.
This applies both to damage caused to the products themselves as a result thereof and to all consequential damage caused thereby.
If the product is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas **FAHL** Medizintechnik-Vertrieb GmbH shall, as far as legally permissible, be free of any liability, including liability for defects.
Should a serious incident occur in connection with this product of Andreas **FAHL** Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.
Sale and delivery of all Andreas **FAHL** Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTCs) which can be obtained directly from Andreas **FAHL** Medizintechnik-Vertrieb GmbH.
The manufacturer reserves the right to make changes to the product at any time without notice.
FAHL® is a trademark and brand of Andreas **FAHL** Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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FAHL® är ett i Tyskland och EU inregistrerat varumärke som tillhör **Andreas FAHL Medizintechnik-Vertrieb GmbH, Köln**.