
FAHL® FOAM

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PICTOGRAM LEGEND

Pictograms listed below you'll find on the product packaging, if applicable.

 REF	Catalogue number		Keep dry
 LOT	Batch code		Do not reuse
	Use by	 MD	Medical devices
	Date of manufacture		Do not use if package is damaged
	Content (in pieces)		CE marking
	Consult instructions for use		Manufacturer
	Keep away from sunlight		

FOREWORD

These instructions for use are valid for **LARYVOX® Foam & SENSOTRACH® Foam**.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Please keep the package for as long as you use the product.

It contains important information on the product!

1. INTENDED PURPOSE

LARYVOX® Foam & SENSOTRACH® Foam are used in combination with a tracheostomy tube or base plates to align the skin level of an irregular tracheostoma that is situated below the skin level.

2. SAFETY INFORMATION

SENSOTRACH® Foam must only be used in combination with a tracheostomy tube.

SENSOTRACH® Foam must never under any circumstances be inserted directly into the tracheostoma!

LARYVOX® Adhesive Foam Disc & LARYVOX® Foam are intended exclusively for use in combination with **LARYVOX® Tape**.

3. INDICATION

The product is indicated for tracheotomised and laryngectomised patients.

4. CONTRAINDICATION

Use of the **LARYVOX® Foam & SENSOTRACH® Foam** should be avoided while undergoing radiotherapy due to the risk of skin irritation.

After undergoing radiotherapy, use **LARYVOX® Foam & SENSOTRACH® Foam** only after consultation with the treating doctor. Start using the tape at the earliest **2 weeks** after completing radiotherapy.

5. COMPLICATIONS AND ADVERSE REACTIONS

The particular adhesive properties of **LARYVOX® Adhesive Foam Disc & SENSOTRACH® Foam PRO** may cause reversible skin irritation in some cases.

The circular opening can under certain circumstances reduce the usable diameter of the tracheostoma, resulting in a reduction of air supply.

6. PRODUCT DESCRIPTION

The speech option is improved by preventing undesired air loss.

The different versions serve for customisation to achieve a uniform skin level.

To accomplish this, the variants can be used in combination.

LARYVOX® Foam

In combination with **LARYVOX® Tape** the contact surface of the tapes is optimised and thus enables longer wearing times.

7. APPLICATION

Clean the skin around the tracheostoma in front of a mirror using warm water and pH-neutral soap or **OPTI-FAHL®** (REF 33210) stoma cleaning wipes.

Dry the skin.

If necessary, remove adhesive residues from the skin with **OPTICLEAR®** (REF 33500) before cleaning.

To protect the skin from irritation, apply **OPTIGARD®** (REF 33600) to the affected area of the skin.

LARYVOX® Foam

After removing the protective foil backing, place the **LARYVOX® Foam** centred over the tracheostoma.

The different sizes are used for individual adjustment to achieve an even skin level for the covering with **LARYVOX® Tape**.

The sizes can be combined for this purpose.

Finally, place **LARYVOX® Tape** over the **LARYVOX® Foam** according to the instructions for use.

Remove the **LARYVOX® Foam** and **LARYVOX® Tape** by carefully removing the adhesive strips from the skin.

CAUTION.

Do not use creams that contain fat! Otherwise, secure attachment of the tape is no longer ensured.

SENSOTRACH® Foam

The tracheostomy tube should be cleaned and, if necessary, disinfected before insertion.

Please observe the manufacturer's instructions with regard to cleaning, disinfection and care of the tracheostomy tube.

SENSOTRACH® Foam PRO

Remove the protective foil on the back. The tracheostomy tube is inserted through the circular opening provided. This is done with the adhesive surface facing toward the tracheostoma.

SENSOTRACH® Foam

The tracheostomy tube is inserted through the circular opening provided.

The different versions (ROUND/OVAL) serve for customisation to achieve a uniform skin level. To accomplish this, the variants can be used in combination.

CAUTION!

Do not use any fatty creams!

8. HYGIENE INSTRUCTIONS

They must not be cleaned or disinfected, as the use of disinfectants can render the **SENSOTRACH® Foam** & **LARYVOX® Foam** unusable.

Sterilisation and reuse can impair the function and is therefore not permissible either.

9. STORAGE

Store in a dry environment away from sunlight and/or heat.

10. SERVICE LIFE

The maximum period of use / wearing period is **24** hours.

If skin irritations occurs, the **SENSOTRACH® Foam/LARYVOX® Foam** must be changed immediately. Deviations from the expected wearing period occur as a result of individual wearing conditions and skin properties.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies both to damage caused to the products themselves as a result thereof and to all consequential damage caused thereby.

When using the product beyond the period of use specified under Section **11**, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTCs) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

LARYVOX® & **SENSOTRACH®** are trademarks and brands of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne, registered in Germany and the European Member States.

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