


REF 58000 TRACHEAL-TASCHENINHALATOR

BILDER

PIKTOGRAMM LEGENDE

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


GRAPHICS

PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

- REF** Catalogue number
- LOT** Batch code
- Use by
- Date of manufacture
- Content (in pieces)
- Consult instructions for use
- Keep away from sunlight
- Keep dry
- Do not reuse
- For reuse on a single patient
- MD** Medical devices
- Non-sterile
- Do not use if package is damaged
- CE** CE marking
- CE 0482** CE marking identification number of the notified body
- Manufacturer

LARYVOX® NEBULIZER

1. FOREWORD

These instructions for use are valid for **LARYVOX® NEBULIZER**.
 These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.
 When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.
 Please read the instructions for use carefully before using the product for the first time!
 Please keep the package for as long as you use the product.
 It contains important information on the product!

2. SAFETY INFORMATION

Do not fill hot liquids into the container.
 Danger of scalding!
 Do not insert the riser tube into the tracheostoma.
 Risk of injury!
 The **LARYVOX® NEBULIZER** is a single-patient product and is thus intended exclusively for one-time or repeated use in one individual patient.

3. INTENDED PURPOSE

The **LARYVOX® NEBULIZER** is a pocket inhaler for generating a spray mist to moisten the mucous membranes of the tracheostoma.
 Suitable for use with water, saline solution as well as medications and other inhalants approved for use in inhalation therapy.

4. INDICATION

The product is indicated for tracheotomized and laryngectomized patients.

5. CONTRAINDICATION

There are no known contraindications for this product.

6. COMPLICATIONS AND ADVERSE REACTIONS

There are no known complications or adverse reactions for this product.

7. APPLICATION

Unscrew the atomiser head of the container.
 Fill the inhalant into the container.
 Screw the atomiser head on the container.
 Hold the riser tube in front of the tracheostoma.
 Pressing the pump atomiser generates a spray mist with which the mucous membranes of the tracheostoma can be moistened.
 For better/more precise positioning of the spray mist, the riser tube on the atomiser head can be adjusted.

8. HYGIENE INSTRUCTIONS

To avoid the formation of bacteria, the inhalant should be changed daily.
 For hygienic reasons we recommend daily cleaning.
 Thoroughly rinse the container, the atomiser head and the riser tube under lukewarm, running water.
 Allow all parts to dry in a clean environment.

CAUTION!

Do not use water with a temperature > 60 °C for cleaning, as this can lead to damage to the product.

9. STORAGE

Store in a dry environment away from sunlight and/or heat.

10. SERVICE LIFE

The service life depends on the frequency of use.
 The products must be replaced immediately when damaged.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.




12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.
 In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.
 This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.
 When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.
 Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.
 Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.
 The manufacturer reserves the right to make changes to the product at any time without notice.
LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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