Instruction for use

TRAKEOSTOMİ HME(İSİ-NEM DEĞİŞTİRİCİ) FİLTRELER TRACHEOSTOMY HME (HEAT-HUMIDITY EXCHANGER) FILTERS

IFU ID: AW.04.ALT.AL.1002-REV.09 (06/2024)


Contents

AR	العربية (تعليمات التشغيل).....	3/32
CZ	Český (Návod K Použití).....	5/32
DA	Dansk (Brugervejledning).....	7/32
DE	Deutsch (Bedienungsanleitung).....	9/32
EL	Ελληνικά (Οδηγισ Χρησης).....	11/32
EN	English (Instructions For Use).....	13/32
ES	Español (Instrucciones De Uso).....	15/32
FR	Français (Mode D'emploi).....	17/32
HU	Magyar (Használati Utasítás).....	19/32
IT	Italiano (Istruzioni Per L'uso).....	21/32
NL	Nederlands (Gebruiksaanwijzing).....	23/32
PL	Polski (Instrukcje Użytkowania).....	25/32
RO	Romana (Instrucțiuni De Utilizare).....	27/32
RU	Русский (Инструкция По Применению).....	29/32
TR	Türkçe (Kullanım Talimatları).....	31/32

Manufactured by:
Meditera Tıbbi Malz. San. ve Tic. A.Ş.
İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
35900 Tire İzmir
TÜRKİYE
Tel: (+90 232) 513 51 10 (Pbx) www.meditera.com.tr

(TRAKEOSTOMİ HME FİLTRESİ) TRACHEOSTOMY HME FILTER

A

- Aspirasyon Portu
- Absaugöffnung
- Port d'aspiration
- Porta di Aspirazione
- Aspiratie poort
- Puerto de Succión
- Sugeport
- Port odsysania
- Sací port
- Leszívó nyílás
- Θυρίδα απορρόφησης
- منفذ المص
- پورت آسپیراسیون
- 吸口
- Portul de Aspirare
- Канал Для Отсасывания

B

- O2 Portu
- Zusätzlicher O2-Anschluss
- Port O2
- Porta O2 supplementare
- O2 Poort
- Supplerende O2 port
- Port dodatkowego tlenu
- Doplnkový port kyslíku O2
- Kiegészítő O2 nyílás
- Θυρίδα O2
- إضافي O2 منفذ
- O2 پورت
- 补充氧气港
- Portul O2
- Переходник для дополнительной подачи кислорода

C

- Hasta Trakeostomi Tüp Bağlantısı
- Verbindung der Patienten- Trachealkanüle
- Connexion du tube de trachéostomie patient
- Connessione della cannula tracheostomica di paziente
- Tracheostomie buis verbinding van patiënt
- Conexión de Tubo de Traqueostomía al Paciente.
- Patient tracheostomi rørforbindelse
- Połączenie z rurką tracheostomijną pacjenta
- Připojení k tracheostomické trubice
- Beteg tracheostomás tubus csatlakozó
- Υποδοχέα προχεοστομίας του Ασθενούς
- وصلة أنبوب مريض القصبة الهوائية
- وصل تیوپ تراکستومی به بی‌مار
- 患者气管切开术管连接
- Conectarea tubului de traheostomie la pacient
- Переходник для соединения с трахеостомической трубкой, установленной пациенту

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Definition:

Tracheostomy HME device is connected to the end of the tracheostomy tube of patients with tracheostomy with spontaneous breathing and is used to humidify the respiratory air with its heat-moisture exchange feature.

Expected Clinical Benefit:

To protect the patients from exposure to dry and cold air by keeping the heat and humidity in their own breathing air and to provide passive humidification.

Artificial Humidifiers (Heat&Moisture Exchanger):

HME (Heat and Moisture Exchanger) filters are used to humidify breathing air. HME devices protect patients from exposure to dry and cold air by keeping the heat and humidity in their own breathing air. Moisture and heat retained during expiration add heat and moisture to the breathing air during inspiration.

They are filters that only humidify with heat and humidity change. It does not have bacteria and virus filtration properties.

This product has been produced in accordance with the requirements of the ISO 9360 standard.

Intended Use:

It is attached to the tracheostomy tube in patients with spontaneous breathing tracheostomy and is used to humidify the breathing air.

Patient Group to Use the Product :

Mechanically ventilated adult and pediatric patients.

Pediatric	50 <Tidal volume < 300 ml.
Adult	Tidal volume>300 ml.

Users :

For this medical device, users are defined as all healthcare professionals including medical doctors, nurses, anaesthesia technicians, paramedics. Defined users have necessary expertise, training and knowledge to use, assemble and install this medical device. This medical device should only be used and assembled by defined users.

WARNING: This product should be applied by only competent health employee.

WARNING: It should not be used in patients with low tidal volume.

Indications:

Tracheostomy HME device is used in patients with tracheostomy. It provides effective humidification for the following purposes.

- To reduce heat loss in bronchopulmonary tree,
- To help prevent changes in lung function and condensations of secretions,
- Supplying additional oxygen when necessary,
- To facilitate the removal of secretions,
- Maintaining humidity and minimizing tracheal epithelial cell damage.

Tracheostomy HME is for use only in adult and pediatric patients over 15 kg.

Contraindications :

Features of HME Filters (Heat&Moisture Exchanger) are contraindicated in patients with fulminant (suddenly occurring and rapidly worsening), foamy secretion in the airways and lungs, and in patients with hemoptysis.

HME (Heat Moist Exchange) filters should not be used together with active humidifying and nebulizers.

Do not add water or steam to HME filters to increase humidification.

Respiratory filters should not be used in patients with very low tidal volume (e.g. premature).

Operation Method:

By making sure the attachment points are secure and secure; connect HME Device to the type of tracheostomy. In case of supplying oxygen, check that the connections of the oxygen hose are also tightly made. Make sure the airway is clear and secure before tying it up.

- When suction is required, open the suction port cover in order to insert the suction catheter.
- Attach the aspiration catheter and begin suctioning.
- When aspiration is complete, remove the aspiration catheter and close the aspiration port cover.

Accessory Use:

Oxygen Tubing:

When used with a filter as an accessory, it is connected between the oxygen source and the supply port. It helps to provide fresh oxygen support to the patient.

Elbow Connector:

When used with a filter as an accessory, it is used between the circuit and the filter, or between the filter and the patient interface, to ensure the appropriate position for the patient.

CAUTION: Make sure to hold the tracheostomy tube in place and stable when exchanging the HME Device.

CAUTION: When the suction port is not used on the tracheostomy HME device, ensure that the aspiration cap, if present, is closed and always check that the port is fully closed to prevent any possible leakage.

CAUTION: If oxygen support is required, the hose of a suitable oxygen supply should be connected to the 5.5mm oxygen port. Oxygen flow should be controlled according to the patient's needs. After connecting the oxygen hose, always check that the hose at the port is completely tight to prevent any possible leakage.

Connect the device securely to breathing system, catheter mount anesthesia/ventilator machine, closed suction system ensuring that the connections are tight and safe. Also check if present perfect seal of the cap for capnometry port.

WARNING: Please perform control tests for any blockage and leakage on the product and all respiratory line before starting the treatment. It should be ensured before usage that the filter is not clogged and air can pass through it.

WARNING: The product is manufactured under clean room conditions. Do not use if the product package is damaged.

WARNING: Additional dead space from the device during therapy must be taken into account.

WARNING: Apply this single-use product as soon as opening the packaging. The user should replace the product regularly in line with hygiene rules. Replace the product according to the hospital infection committee and clinical rules.

WARNING: This product should be used considering dead space, tidal volume, age, weight and the patient's current diseases.

WARNING: Check the interior of the product in terms of liquid or solid secretion. Replace the device without delay in case of a visible liquid or solid accumulation. Otherwise, risk of pressure increases and insufficient patient ventilation occurs.

WARNING: Droplets may be seen in HME (Heat Moist Exchange) filter. This is normal and does not affect filter performance as long as filter integrity is maintained.

WARNING: Avoid contact with chemical substances, cleaning materials or hand sanitizer.

Device Usage Period per Patient :

The filter is produced for use on a single patient. Tracheostomy HME filters should be changed in case of clogging or at most every 24 hours, whichever comes first.

Immediately replace the filters whose resistance has increased as a result of secretions from the patient.

Explanation of Symbols	
	StockCode
	LotNumber
	Quantity
	ExpireDate
	Manufacturer Information
	UniqueDeviceIdentifier
	Single Use
 eIFU indicator	Read the IFU
	Do not use if package is damaged
	Keep away from sunlight
	Caution
	Medical Device
	Temperature Limits
	Humidity Limits
	Keep dry.
	Non sterile.
	Do not contain Latex
	Do not contain phthalate.
	Notified Body Number

WARNING: Apply this single-use product as soon as opening the packaging. The user should replace the product regularly in line with hygiene rules. Replace the product according to the hospital infection committee and clinical rules.

Disposal of Product:

After use, the medical device should be disposed of in accordance with applicable hospital, hygiene and medical waste disposal regulations.

WARNING: Do not reuse this disposable product. Reuse may result in contamination of infectious agents or interruption of treatment, serious harm or death.

WARNING: Do not immerse, wash, clean, disinfect, or sterilize the product. Reprocessing of disposable medical devices may result in performance degradation or loss of function.

Storage Conditions:

Please store in a cool and dry area without exposing of sunlight. Check the containers whether there is any damage. Store the products by considering storage conditions on the labels.

Shelf Life:

If it is stored under appropriate storage conditions, the product may be used until the expiration date indicated on the label.

General Information about the Product:

The product does not contain latex and phthalate.

Shipment of Product:

The products can be placed in a box in a stretched manner according to the agreement with the customer is shipped as a pallet or bulk.

Label Information:

The product label has been prepared with reference to the 2017/745 Medical Device Regulation, EN ISO 15223 and EN ISO 15986 Standard. Explanatory table regarding the symbols on the label is given below.

Explanation of Symbols	
	StockCode
	LotNumber
	Quantity
	ExpireDate
	Manufacturer Information
	UniqueDeviceIdentifier
	Single Use
	Read the IFU
	Do not use if package is damaged
	Keep away from sunlight
	Caution
	Medical Device
	Temperature Limits
	Humidity Limits
	Keep dry.
	Non sterile.
	Do not contain Latex
	Do not contain phthalate.
	Notified Body Number

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