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
ASSKEA GmbH

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99189 Gebesee - Germany

DYSSUCTION® Suction Device M DYSSUCTION® Suction Device S

Instructions for Use

REF 67700 DYSSUCTION® Suction Device S
REF 67900 DYSSUCTION® Suction Device M

Contents

1	Important instructions.....	5
	Instructions in case of transport damage.....	5
	Instructions for customer feedback.....	5
	About these instructions	5
	Notes on the illustrations and symbols in the instructions for use.....	5
2	Intended purpose	6
	Intended Use	6
	Indications	6
	Contraindications.....	6
	Risks and side effects.....	6
	Users	6
	Environment in which the device may be used.....	6
	Environment in which the device must not be used.....	7
3	Safety information	7
4	Guarantee and information on product liability.....	10
	Guarantee.....	10
	Product liability	10
5	Icons and terms	11
5.1	Icons on the device, packaging and accessories.....	11
	Symbols that identify the medical device and accessories	11
	Instructions for working with the medical device and its accessories	11
	Information on environmental, storage and transport conditions	12
5.2	Glossary	12
6	Product description.....	13
6.1	Overview of product and accessories.....	13
	DYSSUCTION® Suction Device M	13
	DYSSUCTION® Suction Device S	13
6.2	Scope of delivery	14
	DYSSUCTION® Suction Device M	14
	DYSSUCTION® Suction Device S	14
6.3	Information on the battery.....	15
6.4	Information on the internal bacteria filter and charcoal filter in the device	16
6.5	Information on the disposable secretion container for the DYSSUCTION® Suction Device M	16
6.6	Information on the disposable secretion container system for the DYSSUCTION® Suction Device S..	17
7	Start-up.....	18
7.1	Start-up of the DYSSUCTION® Suction Device M	18
7.2	Start-up of the DYSSUCTION® Suction Device S	18
7.3	Function test	19

7.4	Connection of the disposable secretion container for the DYSSUCTION® Suction Device M	20
7.5	Connection of the disposable secretion container system for the DYSSUCTION® Suction Device S ..	20
8	Operation.....	23
8.1	Operating and display elements	23
	Operating elements	23
	Display icons	23
8.2	Information on operation and functions.....	23
	Pressure adjustment.....	24
	Language selection	24
	Setting suction and pause times	24
	Patient mode	25
8.3	Operating the suction.....	26
9	Postprocessing.....	28
9.1	After operation with the DYSSUCTION® Suction Device M	28
9.2	After operation with the DYSSUCTION® Suction Device S	28
10	Cleaning and disinfection	30
10.1	WARNINGS.....	30
10.2	Surface of device and holder	30
10.3	Disposable secretion container system DYSSUCTION® Suction Device S	32
11	Maintenance.....	33
12	Troubleshooting and service.....	34
13	Decommissioning and storage	36
	Decommissioning	36
	Storage	36
14	Reprocessing and return	37
	Reprocessing and changing patients in domestic situations	37
	Return for reprocessing	37
15	Disposal.....	37
	Packaging.....	37
	Accessories	37
	Device	37
16	Technical data.....	38
17	Order information	41
	DYSSUCTION® Suction Device M	41
	DYSSUCTION® Suction Device S	41

1 Important instructions

INSTRUCTIONS IN CASE OF TRANSPORT DAMAGE

Please check that deliveries are complete and undamaged immediately on receipt.

All defects must be reported within 3 days, otherwise they cannot be accepted.

INSTRUCTIONS FOR CUSTOMER FEEDBACK

The products described in the instructions for use are being continually developed and improved. For this reason we look forward to receiving all customer feedback, comments and suggestions regarding our products and the accompanying documentation that can contribute to improving the products, our service or the documentation.

ABOUT THESE INSTRUCTIONS

Read the instructions for use in full before using the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** for the first time. Observe the safety instructions in the instructions for use to prevent hazardous situations.

If you require assistance with putting the product into use, operation or maintenance, please contact Andreas Fahl Medizintechnik–Vertrieb GmbH. In addition, please report unexpected operation or incidents to Andreas Fahl Medizintechnik–Vertrieb GmbH.

Keep these instructions where you can consult them at any time.

If you transfer the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** to another person, include the instructions for use with the transfer.

NOTES ON THE ILLUSTRATIONS AND SYMBOLS IN THE INSTRUCTIONS FOR USE

Illustration	Meaning
	Caution – information on risks and how they can be avoided.
	Note with useful information and tips
•	List
1. 2.	Carry out the action in the specified sequence.

2 Intended purpose

INTENDED USE

The **DYSSUCTION® Suction Device M** or **DYSSUCTION® Suction Device S** is a battery-powered, mobile, medical suction device for subglottal suction of secretions with a cuff cannula or a cuff tube with suction attachment.

The **DYSSUCTION® Suction Device M** or **DYSSUCTION® Suction Device S** must never be used on more than one patient at a time!

INDICATIONS

Subglottal suction for patients using a cuff cannula or a cuff tube with suction attachment

CONTRAINDICATIONS

None

RISKS AND SIDE EFFECTS

If the back of the trachea is in contact with the suction opening of the cuff cannula or the cuff tube and causes underpressure over an extended period, this may cause oedema or pressure points.

USERS

Intended users are:

- **medical personnel or professional users,**
who are familiar with and have been trained in how to handle, operate and prepare the device and its accessories and to dispose of the accessories as directed by the instructions for use. Familiarisation and training is provided by Andreas Fahl Medizintechnik–Vertrieb GmbH or an authorised dealer.
- **Relatives of the patient,**
who are familiar with and have been trained in how to handle, operate and prepare the device and its accessories and to dispose of the accessories as directed by the instructions for use. Familiarisation and training is provided by trained and experienced medical personnel or professional users.

ENVIRONMENT IN WHICH THE DEVICE MAY BE USED

- in outpatient and inpatient care (professional healthcare facility environment),
- in the home healthcare environment, particularly for suction of secretion with a cuff cannula or a cuff tube with suction attachment.

ENVIRONMENT IN WHICH THE DEVICE MUST NOT BE USED

- in medical facilities in which equipotential bonding is required (e.g. cardiac surgery).
- in potentially explosive areas.
- in the MRI environment.
- for continuous operation and use outdoors.
- in areas with excessively high humidity and wet rooms.
(e.g. in the vicinity of inhalers or boiling water, bathing or showering).
- in hot areas
(e.g. immediately beside heat sources or exposed to direct sunlight).

3 Safety information**WARNING OF PERSONAL INJURY DUE TO IMPROPER USE**

- For your own protection and the safety of the patient never use the device for other than its intended purpose!
- Do not use the device if damage is visible on the case.
- Do not open the device if you are not trained or authorised for this. Have all services, such as battery replacement, performed by authorised service technicians only.

**WARNING OF PERSONAL INJURY DURING TREATMENT**

- Select the suction setting depending on the patient and the medical indication to avoid potential personal injury.
- We recommend having a spare device available to ensure that treatment can be continued if a device fails.

**WARNING OF PERSONAL INJURY BY STRANGULATION**

Persons may be strangled by tubes or power cords.

- Make sure that unauthorised or other people do not come into the vicinity of the device.
- Remove cables and tubes from the vicinity of the patient when they are no longer required.

**WARNING OF DAMAGE TO HEALTH BY INVOLVEMENT WITH INFECTIOUS OR PATHOGENIC GERMS**

Infectious and pathogenic germs in the suctioned material cause damage to health.

- Always conduct suction with a cuff cannula or a cuff tube with suction attachment.
- Make sure that the suction tube does not come into direct contact with the suction site.
- Clean the device and its accessories as described in these instructions for use.

**WARNING OF NON-APPROVED ACCESSORIES AND SPARE PARTS**

- Use only the accessories and spare parts listed in the instructions for use. The safety and functionality of the device cannot be guaranteed with other accessories.

- Use only clean accessories that are undamaged.

ADDITIONAL PRECAUTIONS FOR SAFE USE IN THE HOME ENVIRONMENT

- Keep the device away from children and pets.
- Make sure that lint and dust are removed from the device.
- It may be difficult to read the display under unfavourable light conditions.

BEWARE OF DAMAGE TO THE DEVICE BY PENETRATION OF LIQUIDS

Protect the device from moisture and clean it only as described in these instructions for use. If the interior of the device comes into direct contact with liquids or solids, it must be inspected by Andreas Fahl Medizintechnik-Vertrieb GmbH.

BEWARE OF DAMAGE TO THE DEVICE BY INCORRECT ELECTRICAL CONNECTION

- Use only the included ATM024T-W120 power supply unit (Art. no. 67912).
- When using the power supply unit, connect the device plug of the power supply unit to the suction device first and only then connect the power supply unit to the mains power.
- The power supply unit must be disconnected from the power supply in exactly the reverse sequence.

BEWARE OF DAMAGE TO THE DEVICE BY IMPROPER HANDLING

- Never suction flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Do not use the device if obvious damage is visible on the case.
- Make sure that the device cannot be subject to impact or fall so persons cannot be hit by the falling device.
- Never touch parts of non-medical devices at the same time in the area of the patient and near the patient.

BEWARE OF DAMAGE TO THE DEVICE BY HEATING

- Do not cover the power supply unit.
- Do not expose the device to direct sunlight.
- Keep the suction device, the power cord and the power supply unit away from other heat sources.

**BEWARE OF INSUFFICIENT DEVICE POWER**

When using the devices at an ambient air pressure below 825 hPa, the maximum vacuum of -300 mbar cannot be reached to 100%.

However, this will not affect the safety and functionality of the device.

Overview of air pressure at 15 °C:

Air pressure	Altitude above mean sea level
1013 hPa	0 m
900 hPa	1000 m
825 hPa	1700 m
746 hPa	2500 m
700 hPa	3000 m

**BEWARE OF RADIO INTERFERENCE**

- Other devices, examinations or treatments may be influenced by the device. For this reason, other devices and examinations or treatments going on at the same time should be carefully monitored to ensure that any interference will be detected as soon as possible.

The device fulfils the requirements of EN 60601-1-2. The system will not be affected by the electromagnetic radiation of other devices that conform to the same standard. The system also does not generate electromagnetic radiation that exceeds the limit values of EN 60601-1-2.

More information on electromagnetic immunity and electromagnetic radiation can be obtained from Andreas Fahl Medizintechnik–Vertrieb GmbH by email to vertrieb@fahl.de or it can be downloaded from our website at www.fahl.com.

4 Guarantee and information on product liability

GUARANTEE

The **DYSSUCTION® Suction Device M** or **DYSSUCTION® Suction Device S** is guaranteed for 2 years. It will not be extended or renewed if repairs under guarantee are required.

The battery is guaranteed for 6 months. Wear parts are not included in the guarantee.

Warrant claims will not accepted if:

- the device is opened by unauthorised persons,
- the security seal is removed or damaged,
- repairs are made by unauthorised persons,
- the device has been modified,

because in these cases the basic safety of the device can no longer be guaranteed and its functions may be restricted.

No liability will be accepted for damage caused by the use of non-recommended accessories and spare parts or non-intended use.

Andreas Fahl Medizintechnik–Vertrieb GmbH guarantees to service the device for at least its service life of 5 years.

PRODUCT LIABILITY

The operator of the device will be liable for the functioning of the device if:

- the device is used for other than its intended purpose,
- the device is not used in accordance with the instructions for use,
- the device is opened by unauthorised persons,
- the security seal is removed or damaged,
- assembly, adjustments, expansions, reprocessing or repairs are conducted by unauthorised persons,
- unapproved accessories and spare parts are used.

5 Icons and terms

5.1 Icons on the device, packaging and accessories

Observe the icons on the device, the packaging and the accessories.

SYMBOLS THAT IDENTIFY THE MEDICAL DEVICE AND ACCESSORIES

Illustration	Meaning	Illustration	Meaning
	Manufacturer		Date of manufacture with country of manufacture
	Distribution partner		The product is a medical device.
	Catalogue number		Serial number
	Batch number		Unique device identification

INSTRUCTIONS FOR WORKING WITH THE MEDICAL DEVICE AND ITS ACCESSORIES

Illustration	Meaning	Illustration	Meaning
	Follow the directions of the instructions for use.		The applied part of the device is classified as Type BF . (Body Floating)
	Observe the instructions for use.		
	Specifies the adjustable operating and rest times.		Protection type: IP33 : The device is protected from the ingress of solid foreign bodies ≥ 2.5 mm diameter and from water spray from all directions, including is the device is tilted 60° from vertical.
	Do not use if packaging is damaged!		
	The product is intended solely for single use.		
	The product is intended for reuse on a single patient.		Power range of the suction device: Low volume flow (< 20 l/min), Average vacuum (between 20 kPa and 60 kPa)
	Power supply unit		

Illustration	Meaning	Illustration	Meaning
	This device must not be disposed of with domestic waste.		This is an electrical device of protection class II. The device is not connected to the protective earth.

INFORMATION ON ENVIRONMENTAL, STORAGE AND TRANSPORT CONDITIONS

Illustration	Meaning	Illustration	Meaning
	Keep away from moisture.		Observe the temperature limits for operation, storage or transport.
	Observe the air pressure limits for operation, storage or transport.		Observe the humidity limits for operation, storage or transport.

5.2 Glossary

A

applied part (Type BF)	Applied parts are classified according to their area of application. Type BF applied parts must be installed disconnected from the earth and are not suitable for direct cardiac applications.
-------------------------------	--

C

Contamination	Contamination means that bacteria and viruses from the suctioned material have come into contact with the device.
----------------------	---

M

ME system	Abbreviation for Medical Electrical system
MRI	Abbreviation for Magnetic Resonance Imaging

O

Oversuctioning	Oversuctioning means that the suctioned material has penetrated the interior of the device.
-----------------------	---

S

Suctioned material	Suctioned material is a general term for secretions, blood, serous liquids or rinse liquids that may occur during suction of the upper respiratory tract.
---------------------------	---

6 Product description

6.1 Overview of product and accessories

DYSSUCTION® SUCTION DEVICE M

Fig. 1

Fig. 2

- A** Disposable secretion container with integrated suction tube
- B** Lock for container
- C** (OK) (On) and (C) (Off) buttons
- D** Display
- E** (Up) and (Down) arrow keys
- F** **DYSSUCTION® Suction Device M** device
- G** Connection for power supply

DYSSUCTION® SUCTION DEVICE S

Fig. 3

- A** "Bag" external container with "OneWay" suction bag (1 l)
- B** "Bag" external container holder
- C** Connecting tube for "OneWay" secretion container system
- D** Display
- E** Control panel (OK) (On), (C) (Off) buttons and (Up) and (Down) arrow keys)
- F** **DYSSUCTION® Suction Device S** device
- G** Connection for power supply

6.2 Scope of delivery

DYSSUCTION® SUCTION DEVICE M

- DYSSUCTION® Suction Device M
- 2 x disposable secretion containers with integrated bacteria filter, charcoal filter, gelling agent and suction tube with stepped connector (art. no. 67910-1)
- the power supply unit and country adapter ATM024T-W120V (art. no. 67912)
- Test protocol according to IEC 62353
- possibly additional components and approved accessories (depending on order, see chapter 17)

DYSSUCTION® SUCTION DEVICE S

- DYSSUCTION® Suction Device S
- Disposable secretion container system comprising (art no. 67714)
 - "Bag" external container (1000 ml)
 - "OneWay" suction bag (1000 ml)
 - "Bag" external container holder
 - Connecting tube for "OneWay" secretion container system
 - Disposable suction tube with stepped connector
 - the power supply unit and country adapter ATM024T-W120V
- Test protocol according to IEC 62353
- possibly additional components and approved accessories (depending on order, see chapter 17)

In order to use the device for its intended purpose, the device must be connected to the suction attachment of the cuff cannula or the cuff tube by the suction tube.

The schematic view below shows the connection of the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** with the patient and accessories.


Fig. 4

6.3 Information on the battery

The **DYSSUCTION® Suction Device M** and the **DYSSUCTION® Suction Device S** devices contain an internal lithium-ion battery, which enables mobile use. The device can be charged and operated with the included power supply unit.

- ✓ Charge the battery completely before the initial use and recharge after the first applications
- ✓ The internal battery is not automatically maintained in an operating status. Charge the device if the battery symbol in the display indicates a low or flat battery.
- ✓ Avoid frequent short charges.
- ✓ Never store the device in a discharged state.
- ✓ After extended storage of the device the battery must be fully recharged again.

The typical service life of the battery is about 300 charging cycles, after which over 80% of the original capacity remains.

The battery is protected from overheating while charging. If the battery temperature is exceeded during charging due to improper ambient conditions, the charging process is temporarily interrupted for cooling. This measure is for safety and protection of the battery.

6.4 Information on the internal bacteria filter and charcoal filter in the device

The internal bacteria filter in the devices is a self-sealing bacteria filter and is effective against bacteria and viruses.

The exhaust compartment of the devices has a filter to neutralise unwanted odour-causing substances from the device exhaust. This filter is a thin fabric coated with activated charcoal. The activated charcoal in the fabric collects the odour-causing particles from the exhaust and neutralises them. This effectively reduces the spread of odours.

Observe the replacement intervals listed in chapter 11.

6.5 Information on the disposable secretion container for the DYSSUCTION® Suction Device M

The disposable secretion container is a container with a connected suction tube. A bacteria filter, a charcoal filter and gelling agent are integrated into the disposable secretion container.

The disposable secretion container with integrated suction tube is a single-use product and must be discarded after use.

BACTERIA FILTER, CHARCOAL FILTER AND OVERFLOW PREVENTION

The hydrophobic bacteria filter integrated into the disposable secretion container is effective against bacteria and viruses. This integrated filter prevents oversuction in the event of an application error or a full container. If liquid reaches the filter, suction is no longer possible and the display shows the "system closed canister full" error message. The suction process is interrupted. The disposable secretion container must be changed.

The activated charcoal filter in the disposable secretion container prevents the spread of odours.

The device has an internal filter in addition to the external filter integrated into the secretion container. The internal filter is a self-sealing bacteria filter and is effective against bacteria and viruses.

GELLING AGENT

Disposable secretion containers that are filled with suctioned material can be safely transported and discarded without leaks by using the gelling agent and closing the tube clamp. Regardless of the suction intervals, the suctioned liquid is thoroughly gelled after an average gel time of 2-5 minutes (depending on the consistency of the suctioned materials).

6.6 Information on the disposable secretion container system for the DYSSUCTION® Suction Device S

The disposable secretion container system consists of the "Bag" external container, the holder for the "Bag" external container, the "OneWay" suction bag, the disposable secretion container system connecting tube and the disposable suction tube with stepped connector.

The "OneWay" suction bag and the disposable suction tube with stepped connector are single-use products and must be disposed of after use.

The "Bag" external container and the disposable secretion container system connecting tube are reusable accessories that are intended for one patient.

In this case observe the replacement intervals listed in chapter 11.

BACTERIA FILTER AND OVERFLOW PREVENTION

A bacteria filter, a charcoal filter and gelling agent are integrated into the suction bag "OneWay". The self-sealing bacteria filter is effective against microorganisms.

If liquid reaches the filter, suction is no longer possible and the display shows the " system closed – canister full" error message. This prevents oversuction. The "OneWay" suction bag must be changed. The activated charcoal filter in the "OneWay" suction bag prevents the spread of odours.

GELLING AGENT

The "OneWay" suction bag that is filled with suctioned material can be safely transported and discarded without leaks by using the gelling agent. Regardless of the suction intervals, the suctioned liquid is thoroughly gelled after an average gel time of 2-5 minutes (depending on the consistency of the secretion).

7 Start-up

7.1 Start-up of the DYSSUCTION® Suction Device M

Check the power supply unit with the power cord and all accessories to ensure that they are not damaged and replace any components that are damaged.

1. Remove the device and accessories from the packaging.
2. Before switching on the device, leave it to acclimatise for about 2 hours at room temperature (approx. 20 °C) to avoid damage to the device.
3. Position the device upright or horizontal (Abb. 5) on a flat, solid surface without tilting it. The device can also be placed beside the patient's bed. Adjustable holders for tubing and rail systems are available for this purpose.


Fig. 5 DYSSUCTION® Suction Device M horizontal

It is the responsibility of the doctor to decide whether the patient's condition enables wearable use.

4. Connect the device plug of the power supply unit to the suction device and only then connect the power supply unit to the mains power.
5. Charge the battery completely before the initial use.
6. Inspect the device and test its functions (chapter 7.3).
7. Connect the disposable secretion container to the device (chapter 7.4).

7.2 Start-up of the DYSSUCTION® Suction Device S

Check the power supply unit with the power cord and all accessories to ensure that they are not damaged and replace any components that are damaged.


1. Remove the device and accessories from the packaging.
2. Before switching on the device, leave it to acclimatise for about 2 hours at room temperature (approx. 20 °C) to avoid damage to the device.
3. Always place the device on a level, secure surface without tilting it.
The device can optionally be attached to tubing and rail systems and also to bed frames.
4. Connect the device plug of the power supply unit to the suction device and only then connect the power supply unit to the mains power.

5. Charge the battery completely before the initial use.
6. Inspect the device and test its functions (chapter 7.3).
7. Connect the disposable secretion container to the device (chapter 7.5).

7.3 Function test

Inspect the device and test its functions without a connected container before using the device for therapy. Proceed as follows to test the functions:

1. Press the OK button for 1 to 2 seconds  to switch on the device. The following start screen appears:

2. After 5 seconds the following image is displayed:

3. Start the pump by pressing the  button and allow the device to run without a container.

An alarm does not appear

"system closed – canister full" alarm

If the "system closed – canister full" alarm appears, the internal filter of the device is blocked and must be replaced by an authorised service technician.

4. Then press a finger on the tube connector and restart the therapy.

"system closed – canister full"
alarm appears after no more
than 5 seconds

An alarm does not appear


If the alarm does not appear after a second test, have the device checked by a service technician.

5. Confirm the alarm with OK.

7.4 Connection of the disposable secretion container for the DYSSUCTION® Suction Device M

Always keep a spare disposable secretion container at hand. This is essential for safe operation!

- A Disposable secretion container (250 ml) with stepped connector
- B Lock for container
- C Suction port
- D **DYSSUCTION® Suction Device M**
- E Guide rail for container

Fig. 6

1. Remove the disposable secretion container (250 ml) with stepped connector (A) from the packaging.
2. Push the container onto the guide rail (E) on the **DYSSUCTION® Suction Device M** until the disposable secretion container audibly clicks onto the lock (B).

7.5 Connection of the disposable secretion container system for the DYSSUCTION® Suction Device S


FUNCTIONAL FAILURE DUE TO COLLAPSING "ONEWAY" SUCTION BAG

A leak in the "Bag" external container or in the cover of the "OneWay" suction bag may allow ambient air to flow into the "Bag" external container. This may cause the "OneWay" external bag to collapse. Make sure,

- that the cover of the "OneWay" suction bag is tightly connected to the "Bag" external container.
- that all connections are tight and correctly connected.
- that the "Bag" external container is not damaged and the T-piece is tight.

Always keep a spare "OneWay" suction bag at hand. This is essential for safe operation!

Names of connections

- A vacuum port
- B patient port

Fig. 7

1. Attach the holder for the "Bag" external container to the device by it pushing into the guide from above.
2. Remove the "OneWay" suction bag from the packaging and extend it completely.


Fig. 8

3. Insert the "OneWay" suction bag into the "Bag" external container. Press the cover down firmly at the edges to ensure that the seal is good.


Fig. 9

4. Attach the pre-assembled connecting tube of the "OneWay" suction bag to the bottom end of the T-piece on the external container.


Fig. 10

5. Now insert the "Bag" external container into the holder. Connect the vacuum port of the device with the corresponding vacuum port of the "Bag" external container (top end of the T-piece). Use the included connecting tube.


Fig. 11

6. Connect the patient port of the "OneWay" suction bag to the suction tube. Make sure that the tube is tightly seated.


Fig. 12

8 Operation

8.1 Operating and display elements

OPERATING ELEMENTS

Operating element	Function
	Up
	Down
	OK (On, Enter)
	Cancel (Off, Back)

DISPLAY ICONS

Illustration	Meaning of the display
	Battery charged
	Battery low
	Battery flat
	The power supply unit is connected.
	Filter runtime reached; the filter must now be replaced by an authorised service technician!
	Suction time
	Pause time
	Alarm OFF The "System closed" alarm is disabled.

8.2 Information on operation and functions

The vacuum is generated by a diaphragm pump. After switching on, the vacuum pump generates a vacuum in the tube system and secretion container, which enables the suctioned material to be removed by suction (through a cuff cannula or a cuff tube with suction attachment). The suctioned material from the patient is moved to the disposable secretion container or disposable secretion container system.

PRESSURE ADJUSTMENT

As soon as the device is switched on the pressure settings can be individually adjusted.

The pressure can be set within a range of -60 mbar to -300 mbar (in 10 mbar steps).


Use pressure settings that have been approved by the treatment doctor.

LANGUAGE SELECTION

German is set as the default language. Proceed as follows to set the language:


1. Press the OK button for 1 to 2 seconds (OK) to switch on the device. The following start screen appears:
2. Press the two arrow keys (▲▼) simultaneously while the start screen is displayed. The menu *Einstellung / Setup* opens.
3. Select the menu *Sprache / Language* with the (▲▼) arrow keys.
4. Confirm your selection with the (OK) button.

SETTING SUCTION AND PAUSE TIMES

The suction and pause times can be set when starting the device. Proceed as follows to set the times:

1. Press the OK button for 1 to 2 seconds (OK) to switch on the device. The following start screen is shown for 5 seconds:
2. Press the (▲▼) arrow keys while the start screen is displayed. The menu *Parameters* opens.
3. Select the Parameter menu with the (▲▼) arrow keys.
4. Confirm your selection with the (OK) button.
5. Now use the (▲▼) arrow keys to set the suction time [sec.]. The minimum suction time is 10 seconds and the maximum 60 seconds (adjustable in 1-second steps, default suction time: 20 seconds).
6. Confirm your selection with the (OK) button.
7. Now use the (▲▼) arrow keys to set the prescribed pause time [min.]. The minimum pause time is 3 minutes and the maximum

60 minutes (adjustable in 1-minute steps, default pause time: 15 minutes).

8. Confirm your selection with the  button.

PATIENT MODE

The patient mode can be selected when starting the device. In the patient mode the patient runtime can be viewed and reset and the "System closed" alarm can be enabled or disabled.

The "System closed" alarm is disabled by default on delivery. It is used to monitor the flow of secretions and to display blockages in the suction tube or the suction attachment. If the alarm is manually activated, the status is saved for the following treatments.

Proceed as follows to select patient mode:

1. Press the OK button for 1 to 2 seconds  to switch on the device. The following start screen is shown for 5 seconds:
2. While the start screen is displayed press and hold the  button and also press the  button. The authorisation screen appears.
3. Use the   arrows keys to enter the "xxxx" code. Press the  arrow key to select the desired number of the code and confirm the input with the  button. Select the additional numbers of the code with the  arrow key and also confirm this input with the  button.



The authorisation code for the patient mode must only be given to specially qualified technicians. Andreas Fahl Medizintechnik-Vertrieb GmbH offers the corresponding training and the authorisation code.


Passwords must be kept confidential to prevent misuse.

4. After authorisation the patient runtime is displayed.
5. Press and hold the  button for 3 seconds to reset the patient runtime to zero.
6. Press the  button to open the activation of the "System closed" alarms in the following screen:
7. Use the   arrow keys to set the checkmark to enable or disable the "System closed" alarm. If the "System closed" alarm is disabled this status is shown on the display as follows:


Overview screen

Suction time

Pause time

8. Exit patient mode by pressing the button.

If the alarm is disabled, blockages in the suction system will not be detected.

Although the pressure during the pause time is also monitored with the alarm disabled and it is automatically pressurised, if a pressure drop that is too low occurs, it will no longer be visible to the user.

8.3 Operating the suction

Danger to patients during application in a domestic environment

The greater mobility of patients in the domestic environment increases the risk the secretions from the subglottal region will not be sufficiently suctioned and as a result the therapy may not be optimal. For this reason comprehensive education and instruction for patients and regular inspection of the suction system by qualified medical personnel is absolutely essential.

For optimal suction of the suctioned material always position the device below the suction site. Make sure that the suction tube does not form a sump and that it is at least at the height of the patient.

Proceed as follows for the suction:

1. Connect the suction tube to the suction device of the cuff cannula or the cuff tube.


The suction tube must never come into direct contact with the suction site.


2. Press the OK button for 1 to 2 seconds to switch on the device. The following start screen is shown for 5 seconds:

3. The following image appears on the display:
(default target: -120 mbar)

4. Set the prescribed pressure value (target value) with the arrow keys.

Always use a pressure setting as low as possible.

5. Start the pump by pressing the  button. Two values are shown on the display.
6. The bar on the top display section moves from left to right and indicates the suction time.
7. The suction time is followed by the pause time. The bar on the top display section now moves from right to left and indicates the suction time.
8. Press the  button again to stop the suction.
9. You are returned to the overview screen:
10. Switch off the device by pressing and holding the  button for 3 seconds.


If you experience a high volume of secretion or you want to suction without a pause during rinsing, press the  button twice at the start of the pause time or as required to skip the pause. Repeat this step if necessary.

9 Postprocessing

Observe the information concerning damage to health by involvement with infectious or pathogenic germs in chapter 3.

9.1 After operation with the DYSSUCTION® Suction Device M

- A Disposable secretion container (250 ml) with stepped connector
- B Lock for container
- C Suction port
- D **DYSSUCTION® Suction Device M**
- E Guide rail

Fig. 13

1. Switch off the device.
2. Close the tube clamp on the suction tube.
3. Disconnect the suction tube from the suction attachment of the cuff cannula or the cuff tube.
4. Press on the lock at the top of the container (B), hold it down and pull the disposable secretion container horizontally from the device.
5. Dispose of the disposable secretion container (250 ml) with stepped connector.

The disposable secretion container (250 ml) with stepped connector is a single-use product and must be discarded after use.

6. Attach a new disposable secretion container to the device as described in chapter 7.4. Make sure that the disposable secretion container is firmly attached to the device.

9.2 After operation with the DYSSUCTION® Suction Device S

1. Switch off the device.
2. Close the tube clamp of the suction tube.
3. Disconnect the suction tube from the suction attachment of the cuff cannula or the cuff tube.
4. Disconnect the connecting tube of the "OneWay" suction bag from the bottom end of the T-piece on the "Bag" external container.


Fig. 14

5. If necessary:
Disconnect the connecting tube from the vacuum port of the "Bag" external container and remove the "Bag" external retainer from the holder.


Fig. 15

6. Remove the "OneWay" suction bag from the reusable "Bag" external container.
7. Dispose of the "OneWay" suction bag with the suction tube.

The "OneWay" suction bag and the suction tube are single-use articles.


Fig. 16

8. Insert a new "OneWay" suction bag.
See chapter 7.5 for information

10 Cleaning and disinfection

10.1 WARNINGS

WARNING OF DAMAGE TO HEALTH BY INVOLVEMENT WITH INFECTIOUS OR PATHOGENIC GERMS

Infectious and pathogenic germs in the suctioned material may cause damage to health.

- As a minimum wear suitable disposable gloves to avoid contact with suctioned material.

Depending on the patient's condition additional protective measures may be required. Observe the requirements of the operator.

- Carry out the following steps carefully in order to achieve the required minimum of microbes.

WARNING OF POSSIBLE INJURY BY ELECTRIC SHOCK

- Switch off the device before cleaning.
- Disconnect the device from the power supply.

CAUTION WITH DAMAGE TO HEALTH WHEN WORKING WITH DISINFECTANTS

- Follow the manufacturer's directions to avoid damage to health.

BEWARE OF DAMAGE TO THE DEVICE BY INCORRECT CLEANING AGENTS

- Cleaning agents and disinfectants may damage the surface of the device.
- Do not use aldehydes or acetone for cleaning and disinfection.
- Frequent cleaning and disinfection may cause slight discolouration of the plastic parts of the case. However, this will not affect the function.

10.2 Surface of device and holder

Clean and disinfect the surfaces of the device and the "Bag" external container holder regularly, but no less than once a week and when changing patients.

RESTRICTIONS WITH PREPARATION	Stop using the device if it shows visible damage or if suctioned material has penetrated the device. Replace the "Bag" external container holder if it shows visible damage or if its functions are affected. There are no other restrictions.
--	---

PREPARATION BEFORE CLEANING	Follow the instructions in chapter 9.2 up to and including section 6. Remove the holder for the "Bag" external container from the device by it pushing upwards along the guide rails. You should now have the device, disconnected from the power supply, in front of you with the holder for the "Bag" external container if installed.									
CLEANING	Clean off visible soiling with a moist but not wet lint-free cloth. Moisten the cloth with plain water (< 40°C). Wipe the complete surface of the device and the holder for the "Bag" external container with light pressure. The holder for the "Bag" external container can be machine-washed and can be heat-sterilised after washing.									
DISINFECTION	Now disinfect the cleaned surfaces. Follow the manufacturer's directions for using the disinfectant. Spray the complete surface of the instrument. Wipe the complete surface. The device and the holder for the "Bag" external container can be disinfected with the alcohol disinfectant ingredient group. Andreas Fahl Medizintechnik–Vertrieb GmbH recommends the following disinfectants: <table><tr><th>SURFACE DISINFECTANTS</th><th>SUBSTANCE GROUP</th><th>TYPE</th></tr><tr><td>Incidin® Liquid (Ecolab)</td><td>Alcohol</td><td>solution ready for use</td></tr><tr><td>BacilloI® AF (Bode Chemie)</td><td>Alcohol</td><td>solution ready for use</td></tr></table>	SURFACE DISINFECTANTS	SUBSTANCE GROUP	TYPE	Incidin® Liquid (Ecolab)	Alcohol	solution ready for use	BacilloI® AF (Bode Chemie)	Alcohol	solution ready for use
SURFACE DISINFECTANTS	SUBSTANCE GROUP	TYPE								
Incidin® Liquid (Ecolab)	Alcohol	solution ready for use								
BacilloI® AF (Bode Chemie)	Alcohol	solution ready for use								
DRYING	Allow the disinfected surface to dry completely.									
MAINTENANCE, INSPECTION AND TESTING	The device should now be free from visible soiling. You can prepare the device for the next application in accordance with the chapter "Start-up".									
STORAGE	If the device is not to be used for an extended period, store it as directed in the instructions in the chapter "Decommissioning and storage".									

10.3 Disposable secretion container system DYSSUCTION® Suction Device S

Clean and disinfect the reusable components of the disposable secretion container system regularly.

- Holder for "Bag" external container: at least once a week or as required
- "Bag" external container: at least once a week or as required
- Connecting tube: at least once a week or as required

The "Bag" external container and the connecting tube must be replaced when changing patients!

The holder for the "Bag" external container can remain with the device when changing patients. The description of the cleaning and disinfection process can be found in chapter 10.2.

RESTRICTIONS WITH PREPARATION	Replace the "Bag" external container and the connecting hose if they show visible damage or if its functions are affected or when changing patients. Replace the "Bag" external container after no more than three months when it is used on one patient. The accessories must be replaced immediately in case of severe soiling with dried or encrusted secretions.
PREPARATION BEFORE CLEANING	Remove the disposable secretion container system and disassemble the components of the system as described in chapter 9.2 up to and including section 6. Dispose of the "OneWay" suction bag.
CLEANING	Clean off visible soiling on the "Bag" external container with a moist but not wet lint-free cloth depending on the degree of soiling. Moisten the cloth with plain water (< 40°C). Clean the surface of the connecting tube with a moist, lint-free cloth. The "Bag" external container can be machine-washed and can be heat-sterilised after washing.

DISINFECTION	You can now disinfect the cleaned surfaces of the accessories. They can be disinfected with the alcohol disinfectant ingredient group. Andreas Fahl Medizintechnik–Vertrieb GmbH recommends the following disinfectants: <table><tr><th>SURFACE DISINFECTANTS</th><th>SUBSTANCE GROUP</th><th>TYPE</th></tr><tr><td>Incidin® Liquid (Ecolab)</td><td>Alcohol</td><td>solution ready for use</td></tr><tr><td>Bacillol® AF (Bode Chemie)</td><td>Alcohol</td><td>solution ready for use</td></tr></table> Follow the manufacturer's directions for using the disinfectant.	SURFACE DISINFECTANTS	SUBSTANCE GROUP	TYPE	Incidin® Liquid (Ecolab)	Alcohol	solution ready for use	Bacillol® AF (Bode Chemie)	Alcohol	solution ready for use
SURFACE DISINFECTANTS	SUBSTANCE GROUP	TYPE								
Incidin® Liquid (Ecolab)	Alcohol	solution ready for use								
Bacillol® AF (Bode Chemie)	Alcohol	solution ready for use								
DRYING	Allow the cleaned or disinfected surfaces to dry completely.									
MAINTENANCE, INSPECTION AND TESTING	The accessory parts should now be free from visible soiling. You can prepare the device for the next application in accordance with chapter 7.									
STORAGE	Store the parts in a clean and dry place ready for the next use.									

11 Maintenance

Do not open the device if you are not trained or authorised for this. Have all services, such as battery or filter replacement, performed by authorised service technicians only.

Keep the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** in a generally clean condition.

Replace reusable accessories that are intended for one patient as soon as

- they show visible damage,
- if the soiling is so bad that they can no longer be cleaned,
- when changing patients.

The internal bacteria filter in the device must be replaced after contact with the suctioned material ("system closed – canister full" display), on expiration of the filter service life ([!]) displayed) or during maintenance or repairs.

The charcoal filter in the device must be replaced during maintenance or repairs, when the internal bacteria filter is changed or every 2 years.

12 Troubleshooting and service

CONTACT OF THE DEVICE INTERIOR WITH LIQUIDS AND SOLIDS

If the interior of the device comes into direct contact with liquids or solids, it must be inspected by Andreas Fahl Medizintechnik–Vertrieb GmbH.




Contact Andreas Fahl Medizintechnik–Vertrieb GmbH if a fault cannot be corrected with the actions described here.

Fault	Possible cause	Remedy
The device cannot be switched on	Battery is flat	Connect power supply unit
	Fault with keypad foil	Please contact your service partner
Suction not possible/no secretion flow	Cannula or suction tube blockage/kinked	• Check cannula position or cuff pressure • Wash or replace cannula or suction tube • Adjust position of device
	Tube clamp closed	Open tube clamp
	Secretion too viscous	Contact your doctor or nursing care personnel
	DYSSUCTION® Suction Device M: Disposable secretion container full	Replace disposable secretion container
	DYSSUCTION® Suction Device S: "OneWay" suction bag (1 l) is full	Replace "OneWay" suction bag (1 l)
	Internal filter is blocked	Please contact your service partner
	Device is still in the <i>Einstellungen</i> mode	Complete selection (see chapter 8.2) start therapy

- The alarms are exclusively technically triggered alarms, because they are identified by the monitoring of device-related variables.
- All alarm messages must be confirmed with OK (except for "internal error").

- High-priority alarm messages show a **red flashing background colour** and the beeper sounds (3x, pause, 2x, 3x, pause, 2x) every 3 seconds.
- Low-priority alarm messages show a **static yellow background colour** and the beeper sounds periodically (2x) every 16 seconds.

Error message	Status	Possible cause	Remedy
	Pump off. Operating mode is terminated	DYSSUCTION® Suction Device M: Disposable secretion container full	Replace disposable secretion container
		DYSSUCTION® Suction Device S: "OneWay" suction bag (1 l) is full	Replace "OneWay" suction bag (1 l)
		If the alarm also appears if the container is not connected, the internal bacteria filter is blocked.	Please contact your service partner
	Pump off. Operating mode is terminated	Battery is flat	Connect power supply unit
	Pump off. Operating mode is terminated	Internal fault	Plug in power supply unit briefly and then disconnect it, if the fault occurs again after a restart after 60 seconds, contact your service partner.
* 	Pump off. Operating mode is terminated	Cannula or tube blocked	• Check cannula position or cuff pressure • Wash or replace cannula or tube
		Secretion flow fault (tube kinked)	Adjust position of device
		Tube clamp closed	Open tube clamp
* The "System closed" alarm only occurs if it was activated in patient mode (see chapter 8.2).			
	Current operating mode continues in the background	Low battery charge status	Connect power supply unit

Error message	Status	Possible cause	Remedy
	"Standby" operating mode (Alarm after 15 minutes)	The therapy was not started	Start therapy
		The device was not switched off	Switch off the device

13 Decommissioning and storage

DECOMMISSIONING

1. Switch off the device after suction by pressing and holding the  button for 3 seconds.
2. Disconnect the power plug from the mains socket and then the power cord or power supply unit from the device.
3. Remove the accessories as described in chapter 9.
4. Clean the accessories that can be used again as described in chapter 10.
5. Clean the surface of the device as described in chapter 10.
6. Store the device in its shipping carton until the next use.

STORAGE

Store the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** under the approved storage conditions:

ambient temperature:	-25 °C to +60 °C
relative humidity:	up to 93%, non-condensing
air pressure:	700 hPa to 1060 hPa

Charge the battery of the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** before storing the device. This ensures that it can be operated at any time.

If the device will not be used for an extended period (approx. 10 months), the battery must be recharged.

14 Reprocessing and return

REPROCESSING AND CHANGING PATIENTS IN DOMESTIC SITUATIONS

Professional hygienic reprocessing by Andreas Fahl Medizintechnik–Vertrieb GmbH is required in the case of changing patients in a domestic environment!

Send the **DYSSUCTION® Suction Device M** or the **DYSSUCTION® Suction Device S** to Andreas Fahl Medizintechnik–Vertrieb GmbH for this purpose.

RETURN FOR REPROCESSING

For your partners and customers, Andreas Fahl Medizintechnik-Vertrieb GmbH offers fast and professional reprocessing and the performance of required tests.

The device must be cleaned before shipping it to Andreas Fahl Medizintechnik–Vertrieb GmbH. See chapter 10! All disposable products must be disposed of. The power supply unit and power cord must be included for testing. Notify Andreas Fahl Medizintechnik–Vertrieb GmbH of the return beforehand.

15 Disposal

PACKAGING

Recycle the device packaging if it is no longer required.

ACCESSORIES

Single-use products must not be reprocessed and reused. Dispose of single-use products correctly.

The following information is for multiple-use products only.

1. Clean and disinfect the multiple-use products.
2. Recycle the disinfected multiple-use products.

DEVICE

Do not dispose of the device with domestic waste.

Clean and disinfect the product.

Dispose of the device at a recycling centre or send it to Andreas Fahl Medizintechnik–Vertrieb GmbH.

16 Technical data

DYSSUCTION® SUCTION DEVICE M REF: 67900-OS

Unit suction capacity* (measured at tube connector)	max. 8 l/min (low flow) * The details may vary depending on the altitude, the atmospheric pressure and the air temperature.	
Vacuum	Medium vacuum -60 mbar to -300 mbar (in 10 mbar increments) Conversion factor: 1 kPa ~ 7.5 mmHg; 1 kPa ~ 10 mbar	
Accuracy of the pressure display	Target pressure > -120 mbar max. Δ 5% Target pressure ≤ -120 mbar max. Δ 10%	
Mains voltage mains unit	In: AC 100 – 240 V~/50 – 60 Hz/580 – 320 mA Out: DC 12 V/2.0 A	
max. load current	2.0 A	
permissible input voltage	12 V	
Power consumption at 12 V	33 W	
Battery, rechargeable	min. 7.2 V; lithium ion battery	
Charging time for flat battery	6 – 7 hours	
Energy of the battery pack	<80 Wh	
Run time in battery operation	approx. 18 h, depending on the load on the motor	
Run time in mains operation	Permanent operation	
Applied part	Type BF	
Protection class in accordance with IEC 60601-1	II	
IP degree of protection	IP33	
Sound emission	Operation: 35 dB (A) High priority alarm: 53 dB (A) Low priority alarm: 52 dB (A)	
Ambient conditions	<i>Transport/storage</i>	
	Ambient temperature: -25 °C to +60 °C Relative humidity: to 93%, non-condensing Air pressure: 700 hPa to 1060 hPa	
	<i>Operation</i>	
	Ambient temperature: +5 °C to +40 °C Relative humidity: 15% to 93%, non-condensing Air pressure: 825 hPa to 1060 hPa	
Device dimensions (H x W x D)	165 mm x 220 mm x 90 mm	

	Relative humidity: Air pressure:	+5 °C to +40 °C 15% to 93%, non-condensing 825 hPa to 1060 hPa
Device dimensions (H x W x D)	290 mm x 259 mm (container : + 100 mm) x 130 mm	
Weight (device)	2.2 kg	
Container volume	1 l	
Suction tube (Article number: 100712)	Disposable suction tube with stepped connector, dia. 4 mm (inside), length: 180 cm, unsterile	

17 Order information

DYSSUCTION® SUCTION DEVICE M

Article number	Article name	PU
67910	Disposable secretion container (250 ml) with stepped connector	30
100705-2	Holder variable DYSSuction Suction Device M	1
67912	Power supply unit	1

DYSSUCTION® SUCTION DEVICE S

Article number	Article name	PU
67710	"Bag" external container (1000 ml)	1
67711	"OneWay" suction bag (1000 ml)	60
67712	"Bag" external container holder	1
67713	Connecting tube for "OneWay" secretion container system	1
67714	Disposable container system (incl. "Bag" and "OneWay")	1
67716	Doppelfiltersystem® filter replacement set	1
67717	Holder variable DYSSuction® Suction Device S	1
100712-1	Disposable suction tube with stepped connector (unsterile)	10

Contacts:

Andreas Fahl

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August-Horch-Str. 4a

51149 Köln - Germany

Phone +49 (0) 22 03/29 80-0

Fax +49 (0) 22 03/29 80-100

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web www.fahl.com

Andreas Fahl

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August-Horch-Str. 4a

51149 Köln - Germany

Phone +49 (0) 22 03/29 80-0

Fax +49 (0) 22 03/29 80-100

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
web www.fahl.com

ASSKEA GmbH

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99189 Gebesee - Germany

DYSSUCTION® Suction Device M DYSSUCTION® Suction Device S

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REF 67700 DYSSUCTION® Suction Device S
REF 67900 DYSSUCTION® Suction Device M

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Lämna apparaten till återvinningscentral för kassering eller skicka den till Andreas Fahl Medizintechnik–Vertrieb GmbH.