FAHL® TRACHEOSPREIZER

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PICTROGRAM LEGEND

	Catalogue number		Do not reuse
	Batch code		Medical devices
	Use by		Sterilized using ethylene oxid
	Date of manufacture		Do not resterilize
	Content (in pieces)		Do not use if package is damaged
	Consult instructions for use		Simple sterile barrier system
	Keep away from sunlight		CE marking identification number of the notified body
	Keep dry		Manufacturer

STERILE DISPOSABLE TRACHEAL DILATOR

FOREWORD

These instructions for use are valid for **FAHL®** Tracheospreizer.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

1. SAFETY INFORMATION

The user must first have been trained by qualified personnel in the safe handling and application.

Application should always take place in front of a mirror.

The **FAHL®** Tracheospreizer must not be inserted too deeply into the tracheostoma.

2. INTENDED PURPOSE

The **FAHL®** Tracheospreizer enables dilation of the tracheostoma.

3. INDICATION

The product is indicated for tracheotomized and laryngectomized patients.

4. CONTRAINDICATION

There are no known contraindications.

5. COMPLICATIONS AND ADVERSE REACTIONS

There are no known complications or side effects.

6. PRODUCT DESCRIPTION

Three spreadable wings on the **FAHL®** Tracheospreizer enable the tracheostoma to be opened evenly during a change of tracheostomy tube.

7. APPLICATION

Check the use-by or expiry date.

Do not use after this date.

Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged.

Do not use the product if the packaging has been damaged.

To open the tracheostoma, insert the **FAHL®** Tracheospreizer slightly into the tracheostoma (see picture 1).

By pressing the handles together, the **FAHL®** Tracheospreizer is opened and the tracheostoma is dilated (see picture 2).

8. HYGIENE INSTRUCTIONS

FAHL® Tracheospreizer are single-patient products and intended for single use only.

They must not be cleaned or disinfected.

9. STORAGE

Store in a dry environment away from sunlight and/or heat.

10. SERVICE LIFE

FAHL® Tracheospreizer are single-patient products and intended for single use only.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section PERIOD OF USE, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

FAHL® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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