

LARYVOX® DUO BRUSH & LARYVOX® STOMA LIGHT

REF 29110



REF 29120



REF 69101 (optional)



- ▶ For gentle and effective daily cleaning of voice prostheses
- ▶ With a brush to clean the channel of the voice prosthesis in the tracheoesophageal fistula and a sponge swab to remove deposits and residual secretions from the tracheal flange of the voice prosthesis
- ▶ **LARYVOX® Duo Brush Extra** with conically shaped brush head - this enables intensive cleaning of the inside of the valve flap and the adjacent oesophageal flange section
- ▶ **LARYVOX® Duo Brush Extra** must be rotated while moving in the voice prosthesis body - otherwise there is a risk of prosthesis dislocation
- ▶ With **LARYVOX® Clip** to allow the fixation of a stoma lamp (**LARYVOX® Stoma Light**, REF 69101 available separately). The combination of brush and lamp in one tool allows effective cleaning with one hand only and frees the other hand for other operations
- ▶ With integrated mirror on the backside of the box to release the patient from carrying extra equipment for daily care
- ▶ **LARYVOX® Stoma Light**, battery driven (battery include in scope of delivery), with bright LED-light to facilitate effective and safe cleaning of the voice prosthesis

SCOPE OF DELIVERY

LARYVOX® Duo Brush REF 29110

6 x voice prosthesis brushes in a case,
Size 06 - 18

When ordering, please state the size after the article number (see table).

or

LARYVOX® Duo Brush Extra REF 29120

6 x voice prosthesis brushes in a case,
one size fits all

Optionally available:

LARYVOX® Stoma light REF 69101

Size 29110-xx	Size voice prosthesis	REF
06	4 - 6	29110-06
08	7 - 8	29110-08
10	9 - 10	29110-10
12	11 - 13	29110-12
18	14 - 18	29110-18

CE 0482