OPTIFAHL®

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PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

	Catalogue number		Do not reuse
	Batch code		Medical devices
	Use by		Do not use if package is damaged
	Content (in pieces)		CE marking
	Consult instructions for use		Manufacturer
	Keep away from sunlight		
	Keep dry		

GENERAL INFORMATION

These instructions are applicable for **OPTICLEAR®**, **OPTIFAHL®** & **OPTIGARD®**.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

1. INTENDED PURPOSE

OPTIFAHL® is a special wipe for gentle and careful external cleaning of dirt particles and excess skin oils from the skin.

OPTICLEAR® is a special wipe for removing remnants of skin adhesives of all kinds and is suitable in particular around the tracheostoma.

OPTIGARD® is a special wipe for preparatory skin care in the use of adhesive plasters/base plates and skin adhesives.

2. SAFETY WARNINGS

Flammable.

Do not bring into contact with open flames.

OPTICLEAR®, **OPTIFAHL®** & **OPTIGARD®** are intended only for external use.

Do not apply on open wounds or mucous membranes.

Avoid contact with the eyes.

Apply only if the skin is intact.

Do not use on small children.

3. INDICATION

The product is indicated for tracheotomised and laryngectomised patients.

4. CONTRAINDICATION

Do not use if the patient is allergic to the material.

5. COMPLICATIONS AND ADVERSE REACTIONS

If skin irritation occurs (e.g. reddening, itching), stop using the product immediately and consult a doctor.

6. PRODUCT DESCRIPTION

OPTICLEAR®:

The non-woven wipe, which can be gripped without slipping, is soaked with an alcohol-base cleaning lotion.

Remnants of adhesives are quickly and gently removed from the skin and absorbed.

OPTICLEAR® is therefore particularly suitable for cleaning the skin after adhesive plasters/base plates have been removed.

OPTIFAHL®:

The special cleaning wipe is soaked with natural plant-based ingredients derived from chamomile, cucumber, sage, fennel, rosemary oil and witch hazel.

For gentle and careful external cleaning of the skin from dirt particles and excess skin oils.

Ideal for preparing the skin when using tracheostoma and stoma articles.

OPTIGARD®:

OPTIGARD® protects against skin irritations, which can be caused by adhesive plasters/base plates or special skin adhesives.

The skin is covered with a thin protective film, which dries after a few seconds.

7. APPLICATION

Check the use-by or expiry date.

Do not use the product after this date.

First clean and degrease the affected part of the skin in the normal way.

OPTICLEAR®

Using gentle pressure, push the cleaning wipe under the edge of the adhesive plaster/base plate, to allow the fluid it contains to penetrate between the adhesive plaster/base plate and the skin.

Repeat the process until the adhesive plaster/base plate has become completely detached harmlessly from the skin.

OPTIFAHL®

Then clean the skin with the **OPTIFAHL®** special wipe and water, then dry it thoroughly.

OPTIGARD®

Then moisten the skin area completely with the **OPTIGARD®** special wipe.

Make sure that the skin is completely dry before applying the adhesive plaster or the base plate or if necessary a skin adhesive.

The film layer does not have to be removed for each change of plaster.

For optimum skin protection, however, we recommend washing off the film layer with a mild cleaning lotion for each change of plaster.

8. STORAGE

Store in a dry environment away from sunlight and/or heat.

9. SERVICE LIFE

OPTICLEAR®, **OPTIFAHL®** & **OPTIGARD®** are single-patient products and are intended for one single use.

Reuse may impair the function and is therefore not permitted.

10. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

11. LEGAL NOTICES

en

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies both to damage caused to the products themselves as a result thereof and to all consequential damage caused thereby.

When using the product beyond the period of use specified under Section 9, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTCs) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

OPTICLEAR®, **OPTIFAHL®** & **OPTIGARD® FAHL®** are trademarks and brands of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne, registered in Germany and the EU member states.

12. INGREDIENTS:

OPTIGARD®

Water, pvp/va copolymer, isopropyl alcohol, peg-40 hydrated ricin oil, perfume.

OPTICLEAR®

Isopropyl alcohol, isopropyl palmitate, perfume

OPTIFAHL®

aqua, peg-35 castor oil, witch hazel, virginiana leaf water, cucumis sativus fruit water, alcohol, propylene glycol, chlorhexidine digluconate, bisabolol, salvia sclarea oil, rosmarinus officinalis oil, foeniculum vulgare dulce fruit oil, phenoxyethanol, ethylhexylglycerin, linalool, farnesol, d-limonene, geraniol, tocopherol.

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