## SECUTRACH® SHOWER GUARD

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# PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

en

|                                                                                   |                              |                                                                                   |                                  |
|-----------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------------------------------------|----------------------------------|
|  | Catalogue number             |  | Keep dry                         |
|  | Batch code                   |  | For reuse on a single patient    |
|  | Use by                       |  | Medical devices                  |
|  | Date of manufacture          |  | Do not use if package is damaged |
|  | Content (in pieces)          |  | CE marking                       |
|  | Consult instructions for use |  | Manufacturer                     |
|  | Keep away from sunlight      |                                                                                   |                                  |

## SECUTRACH® SHOWER GUARD

### 1. FOREWORD

These instructions are valid for the **SECUTRACH®** shower guard.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

### 2. SAFETY INFORMATION

Never hold the water jet directly on the tracheostoma.

The **SECUTRACH®** shower guard must never be used with a full bath and for water sports.

### 3. INTENDED PURPOSE

The **SECUTRACH®** shower guard serves to protect the tracheostoma while showering and washing hair.

### 4. INDICATION

The product is indicated for tracheotomized and laryngectomized patients.

## 5. CONTRAINDICATION

The **SECUTRACH®** shower guard serves exclusively to protect the tracheostoma while showering and washing hair.

It offers no protection while bathing if the tracheostoma is inadvertently immersed under water.

## 6. COMPLICATIONS AND ADVERSE REACTIONS

Entry of water into the tracheostoma, even just a small amount, cause an urge to cough through the natural protective reflexes.

In this event, the user has to leave the shower immediately or immediately interrupt the process of hair washing.

The **SECUTRACH®** shower guard must then be repositioned.

## 7. PRODUCT DESCRIPTION

The **SECUTRACH®** shower guard is tied around the neck and closed with a button fastener.

It is adjustable for all neck sizes.

## 8. APPLICATION

Check the use-by or expiry date.

Do not use after this date.

If a product is damaged, it must be discarded immediately.

Tie the **SECUTRACH®** shower guard around the neck so that the air openings are positioned below the tracheostoma.

The air openings are located in the lower area of the **SECUTRACH®** shower guard rubber cap.

The air flows in through it.

Close the **SECUTRACH®** shower guard with the aid of the button fastener so it seals tightly to the skin.

### CAUTION!

The **SECUTRACH®** largely prevents water entering the tracheostoma.

As a 100% seal cannot be achieved, special caution is advised in use.

## 9. HYGIENE INSTRUCTIONS

The **SECUTRACH®** shower guard must be cleaned thoroughly after each use to prevent functional impairment to the **SECUTRACH®** shower guard and to minimise the risk of infection.

A mild, pH-neutral washing lotion can be used for cleaning.

After cleaning, rinse the **SECUTRACH®** shower guard thoroughly several times with clear water and then dry it with a clean, soft, lint-free cloth

### PLEASE NOTE!

Do not use any sharp-edged cleaning utensils!

Heating to over 65°C, disinfection with boiling water or steam sterilisation are not allowed and lead to damage of the **SECUTRACH®** shower guard.

## 10. STORAGE

Store in a dry environment away from sunlight and/or heat.

## 11. SERVICE LIFE

Do not use the **SECUTRACH®** shower guard if it is damaged or if there are residues/soiling that can no longer be removed by rinsing.

## 12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

## 13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

**SECUTRACH®** is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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