en

FAHL® TRACHEOSTOMY TUBES

INSTRUCTIONS FOR USE

SILVERVENT®
DURAVENT®
DURAVENT® PRO
DURACUFF®
DURACUFF® PRO
TRACHLINE®
LARYNGOTEC®
LARYNGOTEC® PRO
PHONTRACH®
SPIRAFLEX®
FAHL® BIESALSKI
TRACHEOTEC®
TRACHEOSILC®
TRACHEOTEC® SILC
DURATWIX®
DURATWIX® PRO
DURATWIX® JUNIOR
DURAPUR®
DURAVEC PRO

outer cannula

inner cannula

Neck flange

Connectors and adapters

LINGO sieved

SUCTION suction line

CUFF

Connectors and adapters

LINGO fenestrated

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

outer cannula without inner cannula

outer cannula with one inner cannula

outer cannula with two inner cannulas

outer cannula with three inner cannulas

outer cannula including 15 mm-connector (UNI)

outer cannula including 22 mm combi-adapter (KOMBI)

outer cannula including 15 mm-swivel connector (VARIO)

inner cannula with 22 mm-connector (KOMBI)

inner cannula with 15 mm-connector (UNI)

inner cannula with 15 mm-swivel connector (VARIO)

inner cannula with low profile

inner cannula sieved

inner cannula fenestrated

inner cannula flexible

for ventilation

MRT suitable

Catalogue number

Medical device

Keep away from sunlight

For reuse on a single patient

Use by

Content (in pieces)

Sterilized using ethylene oxide

Non-sterile

Cannula lenght adjustable

With suction line (SUCTION)

CUFF

Stoma button

Sieving (LINGO)

Speaking valve (PHON)

Decannulation plug

Obturator (insertion aid)

Tube holder included

Specially for children

Outer cannula fenestrated

Length short

Length medium

Length long

Length XL

Batch code

Consult instructions for use

Keep dry

Do not reuse

Date of manufacture

Do not use if package is damaged

Do not resterilize

Single sterile barrier system

Single sterile barrier system with enclosed protective packaging

Double sterile barrier system

Single sterile barrier system with external protective packaging

CE marking identification number of the notified body

CE marking

Manufacturer

4. CONTRAINDICATION

Elevated intrathoracic pressure.

Never under any circumstances use cannula variants without cuff for patients on mechanical ventilation!

Tracheostomy tubes or base plates with a speaking valve must not be used on any account by laryngectomised (without a larynx) patients, as this could lead to serious complications, even suffocation!

Do not use tracheostomy tubes if the patient is allergic to the material used.

Only use perforated/fenestrated cannula variants during artificial respiration after consultation with the treating doctor/physician.

5. COMPLICATIONS AND ADVERSE REACTIONS

The following complications may occur during use: pressure sores, stenoses, necroses, skin irritations, inflammations, granulation tissue, the urge to cough, bleeding, dysphagia, pulmonary oedema.

If bleeding persists, consult your doctor immediately!

6. PRODUCT DESCRIPTION

Refer to the product label for the complete scope of delivery.

A complete overview of all tracheostomy tubes including technical data and tracheostomy tube accessories can be found at www.fahl.com.

Neck flange

The neck flange or the retainer ring of the tracheostomy tubes has two side eyelets for inserting a tube holder.

Ensure that the FAHL® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

Tracheostomy tubes without side eyelets are placed in a base plate via retainer ring.

CAUTION!

The base tape must be replaced immediately if the adhesive property wears off.

Connectors and adapters

Connectors and adapters are used to connect compatible tube accessories.

The 15 mm standard connector (UNI) enables secure connection with speaking valves and HME.

The 22 mm connector (KOMBI) enables attaching compatible filters and valve systems with a 22 mm adapter.

A rotary connector (VARIO) enables the relief of tensile and rotational forces during ventilation.

Inner cannula

Inner cannulas feature different connectors and can be selected according to requirements.

The inner cannulas can easily be removed from the outer cannula, thus allowing the air supply to be increased quickly if necessary (for instance in case of dyspnoea).

PLEASE NOTE!

The size of the inner cannula must be selected according to the size of the outer cannula.

CAUTION!

Inner cannulas must never be used without outer cannula but must always be affixed to the outer cannula.

WARNING! DURACUFF® & DURACUFF® PRO

Tube versions with an outer and inner cannula may only be used during ventilation with a clip to secure the inner cannula.

Due to the flow dynamics and rheology, tracheostomy tubes without an inner cannula are recommended for ventilated patients.

GENERAL INFORMATION

These instructions for use are valid for FAHL® tracheostomy tubes.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Please keep the package for as long as you use the product.

It contains important information on the product!

1. INTENDED PURPOSE

FAHL® tracheostomy tubes are intended to stabilise the tracheostomy following laryngectomy or tracheotomy.

The tracheostomy tube is designed to keep the tracheostomy open.

Tracheostomy tubes with cuff can be used for mechanical ventilation, including under anaesthesia.

Percutan variants are designed for additional use in percutaneous dilatation tracheostomy.

2. SAFETY INFORMATION

When used for the first time, the products must only be selected needs-oriented, used and inserted by appropriately trained physicians or medical professionals.

Switching to another tracheostomy tube model may only be done at a physician's discretion.

To ensure the best possible fit and breathing situation, a tube adapted to the patient's anatomy should be selected in each case.

After positioning the patient, check that the tube is seated correctly.

Patients must be trained by medical professionals / medical device consultants on how to use the FAHL® tracheostomy tubes safely.

Do not perform any repairs or alterations on a tracheostomy tube or its components.

Use of the product contrary to the described intended use can lead to serious injuries or even death!

We highly recommend keeping at least one replacement tube of the same type and size on hand to ensure a continuous supply.

Tracheostomy tubes are single-patient products and intended for use by one patient only.

FAHL® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation.

Blockage can cause asphyxiation!

When using FAHL® tracheostomy tubes, the products may become infected for instance with yeast (Candida), bacteria or the like, which may lead to the material becoming cloudy and its service life being shortened.

Tracheostomy tubes that contain metal parts must under no circumstances be used during a therapy with exposure to radiation (radiotherapy), as this can lead to serious skin damage and mucosal irritation!

Tracheostomy tubes with metal components must not be used during magnetic resonance imaging scans (MRI) and radiotherapy.

Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices!

There is a fire hazard in connection with oxygen devices and oxygen therapies.

Do not use in the vicinity of sparks or open flames.

Do not grease connections, tubes or other accessories in order to avoid fire hazard.

Decannulation plug

Decannulation may only be carried out under medical supervision.

The use of a decannulation plug in laryngectomised patients and patients with chronic obstructive pulmonary disease is contraindicated.

The decannulation plug must never under any circumstances be used when the tracheostomy tube is in sealed condition!

The decannulation plug may only be used with perforated outer cannula without inner cannula.

3. INDICATION

Tracheostomy tubes without a cuff are indicated in spontaneously breathing patients who are dependent on the use of tracheostomy tubes.

CUFF: Tracheostomy tubes with cuff are indicated after tracheotomies for any reason whenever a seal between the tracheal wall and the tube is required.

XL: Overlong tubes can be indicated particularly in cases of deep-seated tracheal stenoses.

SUCTION: Tracheostomy tubes with a subglottic suction device are used to aspirate the sub-glottic region to minimise the risk of microaspiration.

DYS: This tracheostomy tube is indicated for patients with symptoms of dysphagia.

The tracheostomy tube features a special subglottic suction device through a double hole in the suction tube.

MRT: MRI tracheostomy tubes are indicated for magnetic resonance imaging as well as for radiotherapy.

LINGO: Perforated FAHL® tracheostomy tubes (LINGO) are designed exclusively for tracheotomised patients with a preserved larynx or laryngectomised patients with a shunt valve (voice prostheses) and enable speech.

4. CONTRAINDICATION

Elevated intrathoracic pressure.

Never under any circumstances use cannula variants without cuff for patients on mechanical ventilation!

Tracheostomy tubes or base plates with a speaking valve must not be used on any account by laryngectomised (without a larynx) patients, as this could lead to serious complications, even suffocation!

Do not use tracheostomy tubes if the patient is allergic to the material used.

Only use perforated/fenestrated cannula variants during artificial respiration after consultation with the treating doctor/physician.

5. COMPLICATIONS AND ADVERSE REACTIONS

The following complications may occur during use: pressure sores, stenoses, necroses, skin irritations, inflammations, granulation tissue, the urge to cough, bleeding, dysphagia, pulmonary oedema.

If bleeding persists, consult your doctor immediately!

6. PRODUCT DESCRIPTION

Refer to the product label for the complete scope of delivery.

A complete overview of all tracheostomy tubes including technical data and tracheostomy tube accessories can be found at www.fahl.com.

Neck flange

The neck flange or the retainer ring of the tracheostomy tubes has two side eyelets for inserting a tube holder.

Ensure that the FAHL® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

Tracheostomy tubes without side eyelets are placed in a base plate via retainer ring.

CAUTION!

The base tape must be replaced immediately if the adhesive property wears off.

Connectors and adapters

Connectors and adapters are used to connect compatible tube accessories.

The 15 mm standard connector (UNI) enables secure connection with speaking valves and HME.

The 22 mm connector (KOMBI) enables attaching compatible filters and valve systems with a 22 mm adapter.

A rotary connector (VARIO) enables the relief of tensile and rotational forces during ventilation.

Inner cannula

Inner cannulas feature different connectors and can be selected according to requirements.

The inner cannulas can easily be removed from the outer cannula, thus allowing the air supply to be increased quickly if necessary (for instance in case of dyspnoea).

PLEASE NOTE!

The size of the inner cannula must be selected according to the size of the outer cannula.

CAUTION!

Inner cannulas must never be used without outer cannula but must always be affixed to the outer cannula.

WARNING! DURACUFF® & DURACUFF® PRO

Tube versions with an outer and inner cannula may only be used during ventilation with a clip to secure the inner cannula.

Due to the flow dynamics and rheology, tracheostomy tubes without an inner cannula are recommended for ventilated patients.

Cuff:

In the product versions with cuff, the very thin-walled, high-volume cuff adapts well to the trachea and ensures reliable sealing if inflated correctly.

The cuff is inflated by applying a defined pressure to the cuff via the Luer connection (standardised conical connection) of the inflation tube by means of a cuff pressure gauge.

Suction:

Secretions that have accumulated above the inflated cuff can be removed via the suction opening in the outer tube of the Suction tracheostomy tubes.

The suction tube leading to the outside can be connected to a syringe or a suction device.

PLEASE NOTE!

Subglottic suction using a suction device may only be performed if the suction device is equipped with a vacuum regulator.

CAUTION!

For tracheostomy tubes with a suction option (SUCTION), be careful not to use suction excessively to avoid drying out the subglottic region.

To achieve this, use the lowest suction that can effectively remove secretions.

Continuous suction with a stronger negative pressure than -25 cmH₂O (-25 hPa) can lead to injuries of the trachea.

Lingo / Phon:

Tracheostomy tubes with fenestration or sieving (LINGO) can be used in conjunction with a speaking valve (PHON) for voice formation

CAUTION!

Tracheostomy tubes with a speech option (LINGO) are only recommended for tracheotomised patients with normal swallowing function and laryngectomised patients with a shunt valve (voice prosthesis).

The use of perforated tracheostomy tubes in association with heavy secretion, a tendency to granulation tissue, scabbing or during radiation therapy should be accompanied by regular medical monitoring and compliance with shorter change intervals, as the perforation in the outer tube can promote the formation of granulation tissue and the depositing of secretion residues.

Percutan:

The product variants for percutaneous dilation tracheotomy come with an obturator and a guide catheter.

These tracheostomy tubes are designed for percutaneous dilatation tracheostomy (PDT) and are part of the FAHL® PDT sets.

7. APPLICATION

7.1. PREPARATION

Carefully examine the product to ensure that the packaging and contents have not been tampered with or damaged.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

Damaged tubes must be replaced immediately.

Check the use-by or expiry date.

Do not use after this date.

CAUTION!

In case of an unstable tracheostoma, the airways must always be secured before removing the tracheostomy tube and a replacement tube must be kept ready for insertion.

Cuff:

The cannula and the cuff must be tested for leakages directly before and after every insertion and thereafter at regular intervals (see Fig. 2).

For this purpose, fill the cuff with 20 cmH₂O (19.6 hPa) to 30 cmH₂O (29.4 hPa) and observe whether a spontaneous pressure drop occurs.

Should a sign of leakage be detectable, the tracheostomy tube must under no circumstances be used any more since it is no longer functioning properly.

7.2. Insertion of the tube

Before application, users should clean their hands (see picture 3).

Remove the tube from the package (see Fig. 4).

Before placing the tracheostomy tube, apply the compress and the retaining tape of the tracheostoma tube, if required, to the tracheostomy tube.

To facilitate insertion of the tracheostomy tube, it is advisable to lubricate the outer tube by wiping with a lubricant (see Fig. 6).

During insertion, hold the FAHL® tracheostomy tubes by the neck flange or the retainer ring with one hand (see Fig. 7).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Special aids for dilating the tracheostoma are also available. This allows the tracheostoma to be dilated gently and evenly, for instance also in emergency situations with collapsing tracheostoma (see picture 8).

When using an aid for dilating the tracheostoma, take care to ensure that the tracheostomy tube, and especially the cuff, is not damaged by friction.

Now carefully insert the tube into the tracheostomy during the inspiration phase (while breathing in) while tilting the patient's head slightly back (see Fig. 9).

The neck flange respectively the retainer ring must not be inserted into the tracheostomy.

Straighten your head once the tube has been inserted further into the trachea.

The indwelling time and thus the replacement intervals of the tracheostomy tube must be determined individually in consultation with the treating doctor/physician.

The frequency of changing depends on the patient's individual secretion management.

Inner cannula

CAUTION! SPIRAFLEX®

The mark on the inner cannula must point upwards!

To improve the insertion of the inner cannula into the outer cannula, it should be covered with a medical lubricant (e.g. OPTIFLUID® Stoma Oil Towel REF 331550 or FAHL® OPTIFLUID® Lubricant Gel REF 36105).

To avoid obstruction of the airways, the inner cannula should not remain in the tracheostomy tube for more than one day.

After cleaning or disinfecting and drying the inner cannula, the outer surface of the inner cannula can be made more supple by rubbing with stoma oil or lubricant.

Cuff:

CAUTION!

Always insert the tracheostomy tube with the cuff in fully deflated condition (see picture 5)!

The cuff pressure must be specified by the attending physician based on the individual needs of the patient.

Unless otherwise instructed by the physician, a cuff pressure of 20 cmH₂O (19.6 hPa) to 30 cmH₂O (29.4 hPa) is recommended.

Under no circumstances should the cuff pressure exceed 30 cmH₂O (29.4 hPa).

The correct cuff pressure must be checked regularly, i.e. at least every 6 hours.

If the trachea cannot be sealed with the pressure values specified by the physician, the entire air should be withdrawn again from the cuff and the sealing process should be repeated.

Due to the permeability of the balloon wall for gases, it is normal for the pressure in the cuff to decline slightly over time, but it can on the other hand also rise unintentionally during gas anaesthesia.

PLEASE NOTE!

Always use a cuff pressure gauge or a clean syringe to inflate cuffs, with a subsequent follow-up check of the cuff pressure.

Detach the instruments from the Luer connection of the inflation tube as soon as the cuff has been inflated.

WARNING!

The cuff must never under any circumstances be inflated with excessive amounts of air, since this can lead to damage of the tracheal wall, tears in the cuff with subsequent deflation, or deformation of the cuff, in which case airway obstruction cannot be ruled out.

In patients undergoing artificial respiration, the cuff pressure should not be allowed to drop below the cuff pressure value specified by the doctor/physician in order to prevent unnoticed aspiration.

Unless otherwise specified by the physician, the pressure must not drop below 20 cmH₂O (19.6 hPa).

DURATWIX® JUNIOR:

Unless otherwise instructed by the physician, we recommend a cuff pressure of 20 cmH₂O (19.6 hPa) to 25 cmH₂O (24.5 hPa).

Under no circumstances may the cuff pressure exceed 25 cmH₂O (24.5 hPa).

TRACHEOTEC® SILC:

The inflation volume of the high-pressure cuff must be specified by the physician.

The high-pressure cuff is inflated exclusively using a syringe.

The optimal inflation of the high-pressure cuff has been reached when no leakage can be heard anymore or the ventilator machine no longer indicates a leakage.

PLEASE NOTE!

If the desired sealing is not achieved within the required cuff pressure limits even after repeated attempts, a larger diameter cannula may possibly be indicated.

Lingo / Phon:

To avoid fusion with the mucosa, position cannulae with sieve fenestration at a sufficient distance from the tracheal wall.

7.3. Removing the tube

Connected products, such as speaking valves or HME (heat and moisture exchangers), must first be removed before the FAHL® tracheostomy tube is withdrawn.

WARNING!

If the tracheostomy is unstable, or in emergency situations, the tracheostomy can collapse after withdrawal of the tracheostomy tube, thereby impairing air supply.

A fresh tracheostomy tube must be kept ready for use in such cases and must be quickly inserted if necessary.

A tracheal dilator can be used for temporarily securing the air supply.

The head should be tilted back slightly for removal of the tube.

Grip the tube at the side by the neck flange resp. the retainer ring (see Fig. 11).

Carefully remove the tracheostomy tube with an anchoring motion.

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

Cuff:

The cuff must be emptied before removing the tracheostomy tube (see picture 10).

Never use a cuff pressure gauge to empty the cuff. Always use a syringe for this.

Before deflating the cuff by means of a syringe and withdrawing the tracheostomy tube, the region of the trachea above the cuff must first be cleaned by suctioning off secretions and mucus. It is recommended that the patient be suctioned while at the same time unblocking the tracheostomy tube.

Suctioning is performed by inserting a suction catheter through the cannula tube into the trachea.

Fully deflate the cuff while suctioning at same time.

Proceed very carefully to avoid injury to the mucus membranes.

PLEASE NOTE! SPIRAFLEX®

If the cuff cannot be unblocked, loosen the locking screw of the neck flange.

8. HYGIENE INSTRUCTIONS

DURATWIX® JUNIOR, TRACHEOSILC®, TRACHEOTEC® SILC

Cleaning, disinfection or (re)sterilisation as well as re-use can impair the safety and function of the product and are therefore not permissible.

Obturator:

DURATWIX® JUNIOR, SPIRAFLEX®, TRACHEOSILC®, TRACHLINE®

Insertion aids are single-patient products and are intended for single use only.

They must not be cleaned or disinfected.

DURACUFF®, DURACUFF® PRO, DURACUFF® MOVE, DURAPUR®, DURATWIX®, DURATWIX® PRO, TRACHEOTEC®

The obturator can be cleaned in the same way as the tracheostomy tube.

All other cannulas:

For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the tracheostomy tube at least twice a day, more often during heavier secretion production.

Please note that the personal cleaning schedule, which may also include additional disinfection if necessary, must always be coordinated with the attending physician.

Tracheostomy tubes are considered as belonging to the category of instruments with hollow cavities. Special care must therefore be taken when performing cleaning or disinfection that the tube is fully wetted with the solution used and is unobstructed (see picture 12).

Remove any connected aids before cleaning.

The solutions used for cleaning and disinfection must be freshly prepared every time.

Rinse the cannula cleaning tub as well as the lid and sieve insert thoroughly under running water.

PLEASE NOTE!

Never under any circumstances expose the tubes to boiling water for disinfection.

The tracheostomy tubes must not be cleaned and disinfected using a dishwasher, nor a steam cooker, a microwave oven, a washing machine or similar appliances!

CAUTION!

Cleaning agent and disinfectant residues on the tracheostomy tube can lead to irritations of the mucus membranes or other health impairments.

Inner cannula

The inner cannula must also be removed from the outer cannula.

PLEASE NOTE! SPIRAFLEX®

Do not use cleaning powder or a cleaning brush for the inner cannula!

Cleaning is carried out under running water.

Gentle and careful kneading movements can be used to remove secretion residues.

speaking valve

In tracheostomy tubes with speaking valve, the speaking valve must first be detached from the inner cannula.

The valve itself must not be cleaned with the brush since it could otherwise be damaged or broken off.

PLEASE NOTE!

Silver speaking valves must be removed before cleaning and disinfection and cleaned separately with water.

Cuff:

Do not use a brush to clean the cuff to avoid damaging the cuff!

8.1. Cleaning

A defect of the tracheostomy tube can only be avoided with careful and diligent handling.

PLEASE NOTE!

Only clean the FAHL® tracheostomy tubes with cleaning agents approved by the manufacturer of the tubes.

Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

CAUTION!

Only use cleaning agents if the tube is outside the tracheostoma.

First, rinse the tracheostomy tube thoroughly under running water.

Use only lukewarm water to prepare the cleaning solution.

We recommend using the OPTIBRUSH® CONT cannula cleaning tub (REF 31200) to facilitate the cleaning process.

To do this, hold the sieve insert by the upper edge to prevent contact with and contamination of the cleaning solution.

Only place one tracheostomy tube at a time in the sieve insert of the cannula cleaning tub.

You can place the inner and outer cannula side by side for this procedure.

Cleaning of the outer cannula must only be performed when the cuff is inflated and the safety balloon must be outside the cleaning solution to prevent intrusion of the cleaning solution into the balloon.

Add 1 level scoops of OPTIBRUSH® CLEAN (REF 331110 or REF 331100) cannula cleaning powder or one cleaning bag of the OPTICIT® Clean Set (REF 331180) to the cannula cleaning tub.

Fill the cannula cleaning tub with water (350 ml) and swivel the cannula cleaning tub gently so that the powder begins to dissolve.

Immerse the sieve insert with the cannula in the prepared cleaning solution.

Carefully swivel the cannula cleaning tub back and forth for 5 seconds.

After the immersion time of 15 minutes, rinse the tracheostomy tube thoroughly under lukewarm running water.

No residues of the cleaning agent should remain on the tracheostomy tube.

If necessary, e.g. if stubborn secretions cannot be removed by the cleaning bath, additional cleaning with a special cannula cleaning brush (OPTIBRUSH® (REF 31850) or OPTIBRUSH® PLUS (with fibre attachment) (REF 31855) can prove helpful.

SPIRAFLEX®, SPIRAFLEX® MRT, SPIRAFLEX® XL:

When doing this, the neck flange should be de-arrested and moved to and fro several times to ensure thorough cleaning underneath the locking screw of the neck flange as well

PLEASE NOTE!

Do not use a damaged cleaning brush!

CAUTION!

Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the cannula cleaning brush into the cannula from the cannula lip (see picture 15).

Use the brush as intended and be extremely careful to avoid damage to the tracheostomy tube.

Rinse the tracheostomy tube thoroughly under lukewarm running water.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

FAHL® BIESALSKI, SPIRAFLEX®:

Components such as the adapter, speaking valve, cough cap, decannulation plug, can be cleaned in the same way as the tracheostomy tube.

Suction:

To clean the suction tube of the SUCTION variants, rinse the suction tube several times with clear water using a syringe until no secretion residues are visible.

To remove residual water from the suction tube, inject air through the suction tube with a syringe until no more droplets of water are visible.

8.2. Disinfection

Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).

Please observe the instructions for use for the respective product.

Add 3 level scoops of OPTICIT (REF 331170) cannula cleaning powder or one cleaning bag of the OPTICIT® Clean Set (REF 331180) to the cannula cleaning tub.

Fill the cannula cleaning tub with water (350 ml) and swivel the cannula cleaning tub gently so that the powder begins to dissolve.

Place the cleaned tracheostomy tube (inner and/or outer cannula) in the solution.

Disinfecting the outer cannula with cuff

Tracheostomy tubes with cuff may only be disinfected provided that this is performed and checked with utmost care.

The cuff must always be inflated beforehand.

Carefully swivel the tub back and forth for 5 seconds.

Then leave to stand for 30 minutes.

After the 30 minutes, rinse the tracheostomy tube thoroughly with lukewarm running water.

Allow the tracheostomy tube to dry and store it properly.

PLEASE NOTE!

Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used.

SILVERVENT®:

For disinfection, place the silver cannula in boiling water (100 degrees Celsius or 212 degrees Fahrenheit) for 5 to 10 minutes.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

CAUTION

Allow the tube to cool down before inserting it again to avoid scalds and burns.

8.3. Sterilisation / Autoclaving

Resterilisation is not allowed.

9. STORAGE

Originally packed tracheostomy tubes must be stored in a dry environment away from sunlight and heat.

Cleaned cannulas that are not in use must be stored in a dry environment in a clean container protected from dust, sunlight and heat.

After cleaning tracheostomy tubes, make sure that the cuff has been deflated (unblocked) for storage.

10. SERVICE LIFE

The maximum service life is not extended by using the tube intermittently (e.g. alternately with other tubes in the scope of the cleaning intervals).

If a pathogenic microbe (e.g. MRSA) is present, the period of use/indwelling time is reduced accordingly.

DURATWIX® JUNIOR, TRACHEOSILC®, TRACHEOTEC® SILC:

The maximum period of use must not exceed 29 days.

DURACUFF®, DURACUFF® Pro, DURACUFF® MOVE, DURATWIX®, DURATWIX® Pro, DURATWIX® JUNIOR, FAHL® BIESALSKI, SPIRAFLEX®, TRACHEOTEC®, TRACHLINE®:

The maximum period of use must not exceed 29 days.

The product must be disposed of after 60 cleaning cycles.

DURAPUR®, DURAVEC® PRO, DURAVENT®, DURAVENT® PRO, LARYNGOTEC®, LARYNGOTEC® PRO, PHONTRACH®:

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14

15

When used and looked after correctly, in particular observing the instructions for use, cleaning, disinfection and maintenance provided in the instructions for use, these tracheostomy tubes have an expected service life of three to six months.

The tracheostomy tubes must be replaced at the latest after 6 months of use or 6 months after the packaging has been opened - if damaged or the like, they must be replaced immediately though - the tracheostomy tubes must not be used any longer and must be disposed of correctly. The product must be disposed of after 360 cleaning cycles.

SILVERVENT®:

Tracheostomy tubes made of sterling silver are designed for long-term use.

Send us your tracheostomy tubes for checking and reconditioning on a regular basis, at least every 6 months.

However, the tracheostomy tubes should be replaced at the latest 3 years after first use.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 10, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.



FAHL®, SILVERVENT®, DURAVEC®, DURAVENT®, DURACUFF®, LARYNGOTEC®, PHONTRACH®, SPIRAFLEX®, TRACHEOTEC®, TRACHEOSILC®, DURATWIX®, DURAPUR® and TRACHLINE® are registered trademarks of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne in Germany and the EU member states.

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