


TRACHEOSTOMA SUCTION CATHETER Soft Metric

REF 68000



REF 60700 (optional)



- ▶ Sterile tracheostoma catheters for suctioning tracheal secretions
- ▶ Made of soft, flexible and transparent plastic
- ▶ With central opening and four lateral openings ("eyes") at the tip to avoid adherence to the trachea
- ▶ With grading in 1 cm increments to facilitate control of the insertion depth
- ▶ Wide range of sizes (CH 6-18)
- ▶ Colour coded connectors
- ▶ Available in two different lengths: 25 cm and 50 cm (the short variant is particularly suited for tracheal suctioning)
- ▶ To be used with the TRACHFLOW® Fingertip, to control the suction power. The Fingertip is connected between the suction device and the sterile stoma catheter. One size of the TRACHFLOW® Fingertip fits all sizes of stoma catheters

SCOPE OF DELIVERY

**TRACHEOSTOMA SUCTION CATHETER
Soft Metric S** **REF 68100**

1 x Suction catheter, sterile, atraumatic


or
**TRACHEOSTOMA SUCTION CATHETER
Soft Metric L** **REF 68000**

1 x Suction catheter, sterile, atraumatic

When ordering, please state the size after the item number!

Optional available:
TRACHFLOW® FINGERTIP **REF 60700**
(unsterile)

INFORMATION



CE 0044 STERILE EO (not valid for REF 60700)

Article	Size	Length	Shore	Unit	REF
TRACHEOSTOMA SUCTION CATHETER Soft Metric S	CH 06-18	25 cm	A70	30	68100-xx
TRACHEOSTOMA SUCTION CATHETER Soft Metric L	CH 06-18	50 cm	A70	30	68000-xx