

LARYVOX® TAPE FLEXIBLE


REF 48120



REF 48220



REF 48320



- ▶ Base plate - transparent appearance with water-repellent effect
- ▶ Non-invasive
- ▶ High adhesive strength on the skin around the tracheostoma
- ▶ Suitable for problematic tracheostoma situations
- ▶ Suitable for normal physical activity
- ▶ Suitable for use in combination with shower guard
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (**KOMBI**)
- ▶ Can be combined with special silicone tubes (**LARYNGOTEC® KOMBI CLIP** or **LARYNGOTEC® KOMBI CLIP LINGO**) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval

SCOPE OF DELIVERY

LARYVOX® TAPE FLEXIBLE

15 x Base Plate, individually packed

1 x Instruction for use

round**REF 48120****oval****REF 48220****XL oval****REF 48320**