



FAHL

OPTIBRUSH® / OPTIBRUSH® PLUS / BASIC


REF 31850



REF 31855



REF 31800



- ▶ Special brush for gentle cleaning of tracheostomy tube
- ▶ Bendable and flexible
- ▶ Bristles in dense array around the stainless steel wire base which is twisted to avoid sharp edges which might damage the tracheostomy tube
- ▶ **OPTIBRUSH® (REF 31850)** comes with a rounded tip to avoid scratching of the soft tracheostomy tube material
- ▶ **OPTIBRUSH® PLUS (REF 31855):** The brush tip is surrounded by a soft, wipe-active and absorbent fibre top which further reduces the risk to damage the tube and which ensures a more effective cleaning
- ▶ **OPTIBRUSH® BASIC (REF 31800):** The brush tip is surrounded by a soft, approx. 1 cm wide, wipe-active and absorbent wool tuft which further reduces the risk to damage the tube and which ensures a more effective cleaning
- ▶ Wide size range

SCOPE OF DELIVERY

OPTIBRUSH® REF 31850-xx

4 x OPTIBRUSH® cleaning brushes
1 x Instruction for use

or

OPTIBRUSH® PLUS REF 31855-xx

4 x OPTIBRUSH® PLUS cleaning brushes
1 x Instruction for use

or

OPTIBRUSH® BASIC REF 31800-xx

4 x OPTIBRUSH® BASIC cleaning brushes
1 x Instruction for use

When ordering, please state the size after the item number!

Product	Sizes	Unit	REF
OPTIBRUSH®	04 - 14	4 pieces	31850-xx
OPTIBRUSH® PLUS	05 - 12/14	4 pieces	31855-xx
OPTIBRUSH® BASIC	06/08/10/12	4 pieces	31800-xx