

LARYVOX® TAPE SPECIAL


REF 48145



REF 48245



REF 48345



- ▶ Stable and breathable material; offers the best possible adhesive properties and an unobtrusive appearance
- ▶ High flexibility simplifies handling, especially for the treatment of difficult tracheostomas
- ▶ Non-invasive
- ▶ Adapts comfortably to movements
- ▶ High flexibility and long wearing period
- ▶ Optimal for physical activity
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (**KOMBI**)
- ▶ Can be combined with special silicone tubes (**LARYNGOTEC® KOMBI CLIP** or **LARYNGOTEC® KOMBI CLIP LINGO**) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval

SCOPE OF DELIVERY

LARYVOX® TAPE SPECIAL

15 x Base Plate, individually packed

1 x Instruction for use

round	REF 48145
oval	REF 48245
XL oval	REF 48345