


NECKFIX® ENDO PRO

REF 32920



REF 32921



REF 32931



REF 32932



- ▶ Two-piece fixation set for endotracheal tubes, enables secure and comfortable fixation
- ▶ Centre section: double-sided velour lamination, low-stretch, with pre-glued adhesive strip and slit for secure fixation of the endotracheal tube
- ▶ Neck pad: soft, absorbent and skin-friendly cotton lamination on the inside, foam padding in the centre, velour lamination on the outside for secure fastening of the Velcro fasteners - available in two sizes: 32.5 cm + 47 cm
- ▶ **NECKFIX® ENDO PRO SET:** centre section and neck pad
- ▶ **NECKFIX® ENDO PRO:** Centre section with inserted adhesive strip and 2 Velcro strips
- ▶ **NECKFIX® ENDO PRO Neck pad L**
- ▶ **NECKFIX® ENDO PRO Neck pad XL**
- ▶ Colour: white-turquoise

SCOPE OF DELIVERY

NECKFIX® Endo PRO Set REF 32920

100 x centre section (28 cm) with adhesive strip and 2 Velcro strips, neck pad L, 32.5 cm

or

NECKFIX® Endo PRO REF 32921

8 x 10 centre section (28 cm) with adhesive strip and 2 Velcro strips

or

NECKFIX® Endo PRO Neck Pad L

8 x 10 neck pads L, 32.5 cm

REF 32931

or

NECKFIX® Endo PRO Neck Pad XL

5 x 10 neck pads L, 47 cm

REF 32932