

NEBUFIRST® PRO

REF 50600



- ▶ Powerful compressed air nebulizer for aerosol therapy
- ▶ Aerosol therapy can be used to treat a variety of common respiratory conditions such as asthma, sinusitis, pharyngitis, tonsillitis, laryngitis, bronchitis, bronchiolitis, bronchiectasis, bronchopneumonia and rhinitis.
- ▶ Compact design with rounded shapes
- ▶ For adults and children
- ▶ Robust device for everyday use
- ▶ Easy to clean
- ▶ 2-year guarantee
- ▶ Maintenance-free piston compressor

SCOPE OF DELIVERY


NEBUFIRST® PRO	REF 50600
1 x Inhalation device	
1 x Nosepiece	
1 x Mouthpiece	
1 x Adult mask	
1 x Paediatric mask	
1 x Medicine cup	
1 x Connecting tube	
2 x Replacement air filter	
1 x Instructions for Use	

TECHNICAL DATA

Dimensions	188 x 148 x 92 mm	Working Enviroment	
Weight	1,35 kg	Temperature	10°C – 40°C
Nebulization Rate	≥ 0,25 ml/min	Humidity	30 % RH – 85 % RH
MMAD	approximately 5 µm	Atmospheric Pressure	700 hPa – 1060 hPa
Nebulization Pressure	200 Kpa – 300 Kpa (29 psi – 43,5 psi)	Storage Environment	
Free Flow Rate	4 L/min. – 7 L/min.	Temperature	-20°C - +50°C
Capacity of Medicine Cup	≤ 8 ml	Humidity	30 % RH – 85 % RH
Sound Level	≤ 65 dB	Atmospheric Pressure	500 hPa – 1060 hPa
Power consumption	approximately 150 VA		
Power source	AC 230 V, 50 Hz		
Mode of operation Duty Circle	20 Minuten ON 40 Minuten OFF		
Overvoltage Category	Category II		
Protection against harmful ingress of water or particulate matter	IP21		

CE 0197

NEBUFIRST® PRO



Adult mask



Paediatric mask



Mount with mouthpiece



Mount for nose



Replacement air filter

DEVICE OVERVIEW	REF
① Inhalation device	
② Tube connection	
③ ON/OFF switch	
④ Air filter cassette	
⑤ Medicine cup	
⑥ Air tube	
⑦ Adult mask	50000-05
⑧ Paediatric mask	50000-08
⑨ Mount with mouthpiece	50000-10
⑩ Mount for nose	50000-06
⑪ 2 x Replacement air filter	50000-07



LARYVOX® Mask

FAHL® Nebulizer Set

