

FAHL® PERCUTAN DILATION SET

sterile

REF 19980



- ▶ Set for performing percutaneous dilatation tracheostomy (PDT) according to Ciaglia¹
- ▶ Dilator with hydrophilic coating for good gliding properties
- ▶ Stepless bougienage due to conical shape
- ▶ Haptic structure and extra large handle for safe handling
- ▶ Upper and lower sides slightly flattened for customised pressure distribution
- ▶ Clearly labelled 38 Fr./41 Fr.
- ▶ Guide wire made of kink-resistant Nitinol (Seldinger wire) in dispenser for placement
- ▶ Sterile drape (size 75 x 90 cm)
- ▶ Sterile packed
- ▶ Optionally available:
 - **REF 19990 FAHL® Percutan Seldinger wire:**
Guide wire made of kink-resistant Nitinol
 - **REF 19991 FAHL® Percutan Seldinger wire with guiding catheter:**
Guide wire made of kink-resistant Nitinol with guiding catheter

SCOPE OF DELIVERY

FAHL® PERCUTAN DILATION SET (CIAGLIA)

REF 19980

- 1 x Puncture cannula (14 gauge) with catheter
- 1 x Guide wire with dispenser
- 1 x Predilator
- 1 x Guide catheter with safety lock
- 1 x Dilator with hydrophilic coating
- 1 x Safety scalpel
- 1 x Syringe
- 1 x Sterile drape
- 4 x non-woven compresses
- 1 x Instruction for use

Optional available:

FAHL® Percutan Seldinger wire

REF 19990

FAHL® Percutan Seldinger wire with guiding catheter

REF 19991

CE 0482

¹ Byhahn C, Lischke V, Halbig S et al. Ciaglia Blue Rhino: Ein weiterentwickeltes Verfahren der perkutanen Dilatationstracheotomie. Anaesthesist 2000; 49: 202–206

FAHL® PERCUTAN DILATION SET GRIGGS

sterile

REF 19981



- ▶ Set for performing percutaneous dilatation tracheostomy (PDT) according to Griggs¹
- ▶ Guide wire made of kink-resistant Nitinol (Seldinger wire) in dispenser for placement
- ▶ Stainless steel dilatation forceps with guide channel
- ▶ Secure engagement of the forceps for locking in closed position
- ▶ Sterile drape (size 75 x 90 cm)
- ▶ Sterile packed
- ▶ Optionally available:
 - **REF 19990 FAHL® Percutan Seldinger wire:**
Guide wire made of kink-resistant Nitinol
 - **REF 19991 FAHL® Percutan Seldinger wire with guiding catheter:**
Guide wire made of kink-resistant Nitinol with guiding catheter

SCOPE OF DELIVERY	
FAHL® PERCUTAN DILATION SET GRIGGS	REF 19981
1 x Puncture cannula (14 gauge) with catheter 1 x Guide wire with dispenser 1 x Predilator 1 x Dilator forceps 1 x Safety scalpel 1 x Syringe 1 x Sterile drape 4 x non-woven compresses 1 x Instruction for use	
Optional available:	
FAHL® Percutan Seldinger wire	REF 19990
FAHL® Percutan Seldinger wire with guiding catheter	REF 19991
CE 0482	

REF 19982



- ▶ Set for performing percutaneous dilatation tracheostomy (PDT) according to Ciaglia¹
- ▶ Dilator with hydrophilic coating for good gliding properties
- ▶ Stepless bougienage due to conical shape
- ▶ Haptic structure and extra large handle for safe handling
- ▶ Upper and lower sides slightly flattened for customised pressure distribution
- ▶ Clearly labelled 38 Fr./41 Fr.
- ▶ Guide wire made of kink-resistant Nitinol (Seldinger wire) in dispenser for placement
- ▶ Obturators can be used universally for all cannulas of any type/brand with an inner diameter of 7, 8 or 9 mm
- ▶ Sterile drape (size 75 x 90 cm)
- ▶ Sterile packed
- ▶ Optionally available:
 - **REF 19990 FAHL® Percutan Seldinger wire:**
Guide wire made of kink-resistant Nitinol
 - **REF 19991 FAHL® Percutan Seldinger wire with guiding catheter:**
Guide wire made of kink-resistant Nitinol with guiding catheter

FAHL® PERCUTAN DILATION SET UNIVERSAL (CIAGLIA)

REF 19982

- 1 x Puncture cannula (14 gauge) with catheter
- 1 x Guide wire with dispenser
- 1 x Predilator
- 1 x Guide catheter with safety lock
- 1 x Dilator with hydrophilic coating
- 3 x Obturator (in size 7, 8, 9)
- 1 x Safety scalpel
- 1 x Syringe
- 1 x Sterile drape
- 4 x non-woven compresses
- 1 x Instruction for use

Optional available:

FAHL® Percutan Seldinger wire

REF 19990

FAHL® Percutan Seldinger wire with guiding catheter

REF 19991

C € 0482