

FAHL® Fix Clip

REF 32460



REF 32462



REF 32463



REF 32464



- ▶ Reclosable Velcro strap with metal clip for sorting and securing tubes and cables close to the patient.
- ▶ **FAHL® Fix Clip** contributes to securing and sorting supply and drainage lines as well as monitoring cables. Also suitable within the context of mobilisation and for mobile patients.
- ▶ Enables multiple opening and closing due to the sturdy Velcro fastener
- ▶ Easy repositioning of the supply and drainage lines due to the metal clip
- ▶ Velcro strap suitable for manual labelling, e.g. date of infusion
- ▶ Numerous variants available:
 - **FAHL® Fix Clip 1:** 11 cm long, with single strap
 - **FAHL® Fix Clip 2:** 17.5 cm long, with single strap
 - **FAHL® Fix Clip 2 PRO:** 17.5 cm long, with triple strap
 - **FAHL® Fix Clip 3:** 24 cm long, with single strap

SCOPE OF DELIVERY

FAHL® Fix Clip 1, 2.5 x 11 cm REF 32460
50 x Velcro strap with metal clip

or

FAHL® Fix Clip 2, 2.5 x 17.5 cm REF 32462
50 x Velcro strap with metal clip

or

FAHL® Fix Clip 2 PRO, 2,5 x 17,5 cm REF 32463
50 x Velcro strap with metal clip

or

FAHL® Fix Clip 3, 2.5 x 24 cm REF 32464
50 x Velcro strap with metal clip

