

BLOM-SINGER® ELECTROLARYNX



- ▶ Handy and functional speaking aid for laryngectomised patients
- ▶ Easy to handle due to its low weight
- ▶ The speaking aid is applied to the user's throat and generates sound oscillations, which are transmitted to the pharyngeal, oral and nasal area through a membrane
- ▶ These oscillations are converted through natural speaking motions and form an understandable replacement voice
- ▶ The sound of the Electrolarynx can be customised by adjusting the pitch.
- ▶ The volume control can be individually and optimally adjusted to the user's circumstances via easy-to-operate control knobs
- ▶ The Electrolarynx speaking aid includes a magnetic USB charging cable for easy charging
- ▶ The included carrying strap makes it easier to carry the speaking aid

SCOPE OF DELIVERY

**Blom-Singer® ElectroLarynx
Digital Speech Aid** **REF EL2000**
 1 x Blom-Singer® ElectroLarynx
 1 x Carrying strap
 1 x USB charging cable
 1 x Instructions for Use

Optionally available:
Oral Adapter **REF EL2002**

TECHNICAL DATA

Frequency range	70 - 165 Hz
Volume range	60 - 83 dB
Period of use	up to 30 h with a fully charged battery
Dimensions basic device (L x ø)	123 x 32 mm
Weight	127.6 g
Material	ABS plastic
Battery	3.7 Li-Ion, rechargeable
Capacity	3400 mAh
USB connection	Magnetic micro USB

Andreas Fahl Medizintechnik-Vertrieb GmbH is the exclusive sales partner for Blom-Singer® products in Albania, Austria, Bosnia and Herzegovina, Belgium, Bulgaria, Croatia, Denmark, Estonia, Germany, Hungary, Iceland, Latvia, Lithuania, Luxembourg, Montenegro, the Netherlands, North Macedonia, Norway, Portugal, Romania, Russia, Serbia, Slovenia, Sweden and Switzerland.