

NECKFIX® ENDO LOCK

REF 32940



REF 32940-M28



REF 32945



- ▶ Three-piece fixation set consisting of a fixation tape and two plasters
- ▶ The back of the patient's head and neck remain free
- ▶ Tube holder with double-sided velour lamination, low-stretch, with pre-glued adhesive strip and slit for secure fixation of the endotracheal tube
- ▶ 2 plasters made of Hydrocolloid for placement on the patient's face
- ▶ Simple Velcro fastening of the tape to the plaster
- ▶ No protruding Velcro edges on the plaster
- ▶ Colour: white-turquoise

SCOPE OF DELIVERY

NECKFIX® Endo Lock-Set **REF 32940**

8 x 10 Pad (28 cm) with adhesive strip and 2 fixing plasters

or

NECKFIX® Endo Lock **REF 32940-M28**

8 x 10 Pad (28 cm) with adhesive strip

or

NECKFIX® Endo Lock Fixing Plasters

80 x Fixing plasters **REF 32945**