

FAHL® TRACHEOSPREIZER

sterile

REF 33502



- ▶ The Tracheospreizer is used for the dilation of the tracheostoma of tracheostomized and laryngectomized patients
- ▶ Also to be used during cannula changing and for emergency intubation, e.g. decannulation
- ▶ The special shape of the rounded tips minimizes the risk of injury during insertion
- ▶ The three spreadable wings allow the opening of the dilator in all directions and thereby ensure an even circular dilation of the tracheostoma
- ▶ Ergonomic design and handling: The dilator opens when the handles are pressed together
- ▶ This is conform to the natural ergonomics of manual handling and provides a more sophisticated and fine-tuned control of the dilation process

SCOPE OF DELIVERY

FAHL® Tracheospreizer
1 x tracheal dilator

REF 33502

CE 0482