

TRACHFLOW® SUCTION CATHETER Soft Metric

REF 68730



REF 68630

REF 60700
(optional)

- ▶ Atraumatic tracheostoma catheters for particularly gentle suctioning of tracheal secretions
- ▶ With rounded tip to minimize the risk of tracheal injury
- ▶ Made of soft, flexible and transparent plastic
- ▶ With central opening and four lateral openings ("eyes") at the tip to avoid adherence to the trachea
- ▶ With grading in 1 cm increments to facilitate control of the insertion depth
- ▶ Wide range of sizes (CH 08 -16)
- ▶ Colour coded connectors
- ▶ Length: 30 cm and 50 cm
- ▶ The short variant is particularly suited for tracheal suctioning
- ▶ To be used with the TRACHFLOW® Finger-Tip, to control the suction power
- ▶ The TRACHFLOW® Finger-Tip is connected between the suction device and the sterile stoma catheter.
- ▶ One size fits all sizes of stoma catheters.

SCOPE OF DELIVERY

**TRACHFLOW® SUCTION CATHETER
SOFT METRIC S** **REF 68630**

- ▶ 100 x Suction catheter, sterile, atraumatic

or

**TRACHFLOW® SUCTION CATHETER
SOFT METRIC L** **REF 68730**


- ▶ 100 x Suction catheter, sterile, atraumatic

When ordering, please state the size after the item number!

Optionally available:

TRACHFLOW® FINGERTIP **REF 60700**
(unsterile)

INFORMATION



08 10 12 14 16

CE 0123 STERILE EO (not valid for REF 60700)

Product	Size	Length	Shore	Color	Unit	REF
TRACHFLOW® SUCTION CATHETER Soft Metric S	CH 10	30 cm	A65		100	68630-10
TRACHFLOW® SUCTION CATHETER Soft Metric S	CH 12	30 cm	A65		100	68630-12
TRACHFLOW® SUCTION CATHETER Soft Metric S	CH 14	30 cm	A65		100	68630-14
TRACHFLOW® SUCTION CATHETER Soft Metric S	CH 16	30 cm	A65		100	68630-16
TRACHFLOW® SUCTION CATHETER Soft Metric L	CH 08	50 cm	A65		100	68730-08
TRACHFLOW® SUCTION CATHETER Soft Metric L	CH 10	50 cm	A65		100	68730-10
TRACHFLOW® SUCTION CATHETER Soft Metric L	CH 12	50 cm	A65		100	68730-12
TRACHFLOW® SUCTION CATHETER Soft Metric L	CH 14	50 cm	A65		100	68730-14
TRACHFLOW® SUCTION CATHETER Soft Metric L	CH 16	50 cm	A65		100	68730-16