

# LARYVOX® DUO BRUSH & LARYVOX® STOMA LIGHT

REF 29110



REF 29120



REF 69101 (optional)



- ▶ For gentle and thorough daily cleaning of voice prostheses
- ▶ Includes a brush for cleaning the inner channel of the voice prosthesis in the tracheoesophageal fistula and a foam sponge for removing debris and secretions from the tracheal flange of the voice prosthesis
- ▶ **LARYVOX® Duo Brush Extra** with conical brush head – it ensures intensive cleaning of the inside of the valve flap and the adjacent oesophageal section of the flange
- ▶ Move the **LARYVOX® Duo Brush Extra** inside the body of voice prosthesis only with rotary movements - otherwise there will be a danger of dislocating the prosthesis
- ▶ Includes the **LARYVOX® Clip** for attaching a stoma light (**LARYVOX® Stoma Light**, available separately). The combination of brush and lamp allows the cleaning process to be performed with one hand while the other hand remains free
- ▶ The mirror on the back of the box enables use by patients themselves for convenient cleaning. It is no longer necessary to keep a separate mirror at hand
- ▶ **LARYVOX® Stoma Light**, battery-powered (battery included), with very bright LED light for effective and safe cleaning of the voice prosthesis

## SCOPE OF DELIVERY

**LARYVOX® Duo Brush**      **REF 29110-xx**  
6 x voice prosthesis brushes in the pouch, sizes 06 - 18

When ordering, please state the size after the item number (see table).

or  
**LARYVOX® Duo Brush Extra**      **REF 29120**  
6 x voice prosthesis brushes in the pouch, single size

Optionally available:  
**LARYVOX® Stoma light**      **REF 69101**

| Size<br>29110-xx | Voice prostheses size |               | REF             |
|------------------|-----------------------|---------------|-----------------|
|                  | Length<br>[Fr]        | Width<br>[mm] |                 |
| <b>06</b>        | 4 - 6                 | 16-24         | <b>29110-06</b> |
| <b>08</b>        | 7 - 8                 | 16-24         | <b>29110-08</b> |
| <b>10</b>        | 9 - 10                | 16-24         | <b>29110-10</b> |
| <b>12</b>        | 11 - 13               | 16-24         | <b>29110-12</b> |
| <b>18</b>        | 14 - 18               | 16-24         | <b>29110-18</b> |
| <b>-</b>         | 4-10                  | 20-24         | <b>29120</b>    |

CE 0482