


LARYVOX® HME High Flow

REF 49810



- ▶ Special wide-pore polyurethane foam to reduce re-spiration resistance
- ▶ Can hence be used up to 24 hours a day, even during physical activities
- ▶ Skin colour to provide inconspicuous design
- ▶ Lightweight; low profile
- ▶ Stable housing – can be closed by light finger occlusion to enable laryngectomized patients
- ▶ With a voice prosthesis to speak
- ▶ With 22 mm-connector (**KOMBI**)

REF 49803



SCOPE OF DELIVERY

LARYVOX® HME HighFlow REF 49810

30 x Filter cassette, individually packed
1 x instruction for use

LARYVOX® O₂ HME HighFlow REF 49803

30 x Filter cassette, individually packed
1 x instruction for use