


LARYVOX® TAPE HYDROSOFT

REF 48130



REF 48230



REF 48330



- ▶ Base plate - particularly skin-friendly
- ▶ Non-invasive
- ▶ High fixation comfort
- ▶ Ability to enclose fluids/secretions
- ▶ Not suitable for use in combination with shower guard!
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (**KOMBI**)
- ▶ Can be combined with special silicone tubes (**LARYNGOTEC® KOMBI CLIP** or **LARYNGOTEC® KOMBI CLIP LINGO**) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval

SCOPE OF DELIVERY

LARYVOX® TAPE HYDROSOFT

15 x Base Plate, individually packed

1 x Instruction for use

round**REF 48130****oval****REF 48230****XL oval****REF 48330**