

LARYNGOTEC® KOMBI CLIP/ LARYNGOTEC® KOMBI CLIP LINGO

REF 114830



REF 114820



- ▶ Tracheostomy tube for long-term use, can be used up to 3-6 months
- ▶ Medical grade silicone: soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ With 22 mm connector (**KOMBI**) to connect HMEs
- ▶ Bulging ring (**CLIP**) behind the 22 mm connector allows connection to adhesive base plates
- ▶ Each size in different lengths
- ▶ With sieve-fenestration **LARYNGOTEC® KOMBI CLIP LINGO (REF 114820)** or without sieve-fenestration **LARYNGOTEC® KOMBI CLIP (REF 114830)**
- ▶ Custom-made sieve-fenestration possible

SCOPE OF DELIVERY

LARYNGOTEC® KOMBI CLIP REF 114830

1 x Tracheostomy tube
1 x LARYVOX® Tape STANDARD round
1 x product ID card
1 x Instruction for use

or

LARYNGOTEC® KOMBI CLIP LINGO REF 114820

1 x Tracheostomy tube with sieve-fenestration
1 x LARYVOX® Tape STANDARD round
1 x product ID card
1 x Instruction for use

When ordering, please state the size and length (OL) after the item number!



KOMBI



LINGO




MRT

CE 0482


LARYNGOTEC® KOMBI CLIP/ LARYNGOTEC® KOMBI CLIP LINGO

Technical Data

REF 114830



REF 114820

**Legend**

All lengths in mm
 OD = Outer diameter
 ID = Inner diameter
 OL = Length over the outer curve
 NL = Length over the middle curve
 L1 = Start of phonation
 L2 = Length of phonation
 BA = Bending Angle

LARYNGOTEC® KOMBI CLIP

REF 114830-xxxx

Size	OD	ID	NL	OL	BA	REF
07	10,2	7,0	38,0	44,0	66,0	114830-0744
			38,0	44,0	59,0	114830-0844
			54,0	62,0	82,0	114830-0862
09	12,4	9,0	36,0	44,0	56,0	114830-0944
			55,0	62,0	77,0	114830-0962
			65,0	75,0	93,0	114830-0975
10	13,4	10,0	36,0	44,0	49,0	114830-1044
			55,0	62,0	67,0	114830-1062
			75,0	85,0	90,0	114830-1085
11	14,4	11,0	36,0	44,0	48,0	114830-1144
			55,0	62,0	66,0	114830-1162
			76,0	87,0	91,0	114830-1187
12	15,4	12,0	36,0	44,0	48,0	114830-1244
			55,0	62,0	66,0	114830-1262
13	16,6	13,0	37,0	44,0	47,0	114830-1344
			55,0	62,0	65,0	114830-1362

LARYNGOTEC® KOMBI CLIP LINGO

REF 114820-xxxx

Size	OD	ID	NL	OL	BA	L1	L2	REF
07	10,2	7,0	38,0	44,0	66,0	12,0	15,0	114820-0744
			54,0	62,0	92,0	16,0	20,0	114820-0762
08	11,2	8,0	38,0	44,0	59,0	12,0	15,0	114820-0844
			54,0	62,0	82,0	16,0	20,0	114820-0862
09	12,4	9,0	36,0	44,0	56,0	12,0	20,0	114820-0944
			55,0	62,0	77,0	16,0	26,0	114820-0962
10	13,4	10,0	36,0	44,0	49,0	12,0	20,0	114820-1044
			55,0	62,0	67,0	16,0	26,0	114820-1062
			75,0	85,0	90,0	16,0	26,0	114820-1085
11	14,4	11,0	36,0	44,0	48,0	12,0	20,0	114820-1144
			55,0	62,0	66,0	16,0	26,0	114820-1162
12	15,4	12,0	36,0	44,0	48,0	12,0	20,0	114820-1244
			55,0	62,0	66,0	16,0	26,0	114820-1262
13	16,6	13,0	37,0	44,0	47,0	12,0	20,0	114820-1344
			55,0	62,0	65,0	16,0	26,0	114820-1362

PIKTOGRAMME / PICTOGRAMS / PICTOGRAMMES / PICTOGRAMMEN / PICTOGRAMAS / BILDER



UNI

de	15 mm-Konnektor (UNI)
en	15 mm-connector (UNI)
fr	Connecteur 15 mm (UNI)
nl	15 mm-connector (UNI)
es	Conector de 15 mm (UNI)
sv	15 mm anslutning (UNI)



VARIO

de	15 mm-Drehkonnektor (VARIO)
en	15 mm-swivel connector (VARIO)
fr	Connecteur rotatif 15 mm (VARIO)
nl	15 mm-swivel-connector (VARIO)
es	Conector rotatorio de 15 mm (VARIO)
sv	Ytterkanyl inklusive 15 mm vridanslutning (VARIO)



KOMBI

de	22 mm-Konnektor (KOMBI)
en	22 mm-connector (KOMBI)
fr	Connecteur 22 mm (KOMBI)
nl	22 mm-connector (KOMBI)
es	Conector de 22 mm (KOMBI)
sv	22 mm anslutning (KOMBI)



LINGO

de	Fensterung oder Siebung (LINGO)
en	Fenestration or sieving (LINGO)
fr	Fenêtrage ou perforations (LINGO)
nl	Vensteren of zeven (LINGO)
es	Fenestration o perforación (LINGO)
sv	Fönsteröppning eller siktning (LINGO)




FLAT-PHON

de	Flachprofil mit Sprechventil (FLAT-PHON)
en	Flat profile with speaking valve (FLAT-PHON)
fr	Profilé plat avec valve de phonation (FLAT-PHON)
nl	Vlak profiel met spreekventiel (FLAT-PHON)
es	Perfil plano con válvula fonatoria (FLAT-PHON)
sv	Platt profil med talventil (FLAT-PHON)



PHON

de	Sprechventil (PHON)
en	Speaking valve (PHON)
fr	Valve de phonation (PHON)
nl	Spreekventiel (PHON)
es	Válvula fonatoria (PHON)
sv	Talventil (PHON)



CUFF

de	Mit CUFF
en	With CUFF
fr	Avec ballonnet (CUFF)
nl	Met CUFF
es	Con balón (CUFF)
sv	Med CUFF



SUCTION

de	Mit Absaugvorrichtung (SUCTION)
en	With suction line (SUCTION)
fr	Avec dispositif d'aspiration (SUCTION)
nl	Met afzuigkanaal (SUCTION) voor subglottisch afzuigen van secret boven de cuff
es	Con línea de succión (SUCTION)
sv	Med sugslang (SUCTION)

O₂

de	Mit einem O ₂ -Anschluss
en	With oxygen support
fr	Avec raccord O ₂
nl	Met zuurstofaansluiting (O ₂ -support)
es	Con puerto de oxígeno
sv	Med O ₂ -anslutning



MRT

de	MRT geeignet
en	MRT suitable
fr	IRM-compatible
nl	Geschikt voor MRT
es	Compatible con resonancia Magnética
sv	MR-kompatibel

de Falls zutreffend, sind diese Piktogramme auf dem Datenblatt zu finden.
 en If applicable, these pictograms can be found on the data sheet.
 fr Le cas échéant, ces pictogrammes se trouvent sur la fiche technique.

nl Indien van toepassing zijn deze pictogrammen te vinden op het gegevensblad.
 es Si procede, estos pictogramas se encuentran en la ficha técnica.
 sv I förekommande fall finns dessa piktogram på databladet.