

HUMIDOPHONE® PLUS RED / HUMIDOPHONE® PLUS O₂ RED

REF 46487-03



CE

REF 46489-03



CE 0482

- ▶ Speaking valve with integrated filter
- ▶ The red colour distinguishes it unmistakably from other attachments, such as a HME
- ▶ The membrane is closed in the basic position (valve close position) and opens during inhalation.
- ▶ Advantage: a high blowing pressure is not required to close the valve to be able to speak. Less fatigue for the patient, better speaking quality.
- ▶ Flat profile, low weight
- ▶ Open-pore polyurethane filter
- ▶ Filter in proximal position (toward the patient)
- ▶ With multi-adapter: fits a 15 mm connector (UNI) and also a 22 mm combi-adapter (KOMBI)
- ▶ **HUMIDOPHONE® PLUS O₂ RED** (REF 46489-03) with oxygen connection (5 mm plug connector)

SCOPE OF DELIVERY

- ▶ Instruction for use

HUMIDOPHONE® PLUS RED

REF 46487-03

30 x speaking valve

or

HUMIDOPHONE® PLUS O₂ RED

REF 46489-03

30 x speaking valve

Product	Unit	REF
HUMIDOPHONE® PLUS RED	30 pieces, in 5 piece packs	46487-03
HUMIDOPHONE® PLUS O₂ RED	30 pieces, in 5 piece packs	46489-03

The HUMIDOPHONE® speaking valve must be only used by patients with a partially or fully preserved larynx!