

LARYVOX® TAPE SPECIAL PLUS

REF 48347



- ▶ Stable and breathable material; offers the best possible adhesive properties and an unobtrusive appearance
- ▶ High flexibility simplifies handling, especially for the treatment of difficult tracheostomas
- ▶ Non-invasive
- ▶ Adapts comfortably to movements
- ▶ High flexibility and long wearing period
- ▶ Optimal for physical activity
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (**KOMBI**)
- ▶ Can be combined with special silicone tubes (**LARYNGOTEC® KOMBI CLIP** or **LARYNGOTEC® KOMBI CLIP LINGO**) if additional tracheal stabilisation is needed

SCOPE OF DELIVERY

LARYVOX® TAPE SPECIAL PLUS

15 x Base Plate, individually packed

1 x Instruction for use

REF 48347

