


TRACHEOTEC® PRO VARIO CUFF SUCTION

2 IC

REF 18932



- ▶ Tracheostomy tube with two inner cannulas (IC), with low-pressure **CUFF** and integrated **SUCTION** line to remove secretions above the cuff
- ▶ Soft, flexible cannula tag with imprint of the essential specifications of the cannula
- ▶ With 15 mm swivel connector (**VARIO**) to buffer tensile or torsional forces that may be generated by the respiratory tube
- ▶ With X-ray contrast strip
- ▶ Available in sizes 6.5 / 7 / 8 / 9 and 10
- ▶ Set includes tube holder, an obturator and cannula ID card
- ▶ Inner cannulas (**REF 18933**) can also be ordered separately, size must be selected according to the **TRACHEOTEC® PRO** outer cannula

SCOPE OF DELIVERY

TRACHEOTEC® PRO
VARIO CUFF SUCTION 2 IC **REF 18932-xx**
 1 x Tracheostomy tube with CUFF and SUCTION line
 2 x Inner cannula (IC)
 1 x Tube holder
 1 x Obturator
 1 x Cannula ID card
 1 x Instructions for Use

Please state the size after the article number when ordering! You will find the size chart in the appendix.



VARIO



CUFF



SUCTION


CE 0482

TRACHEOTEC® PRO VARIO CUFF SUCTION

2 IC

TECHNICAL SKETCH

REF 18932



TRACHEOTEC® PRO VARIO CUFF SUCTION

REF 18932-xx

Size	OD1	ID1	CD	NL	OL
065	9,3	6,5	22,0	65,8	73,5
07	10,0	7,0	30,0	68,7	77,0
08	11,3	8,0	30,0	78,5	87,9
09	12,7	9,0	30,0	84,4	94,9
10	14,0	10,0	30,0	93,7	105,3

Legend

All specifications in mm

OD1 = Outer diameter at the neck flange

ID1 = Inner diameter at the cannula tip of the outer cannula

CD = CUFF diameter

NL = Length over the central arch

OL = Length over the outer curve

PIKTOGRAMME / PICTOGRAMS / PICTOGRAMMES / PICTOGRAMMEN / PICTOGRAMAS / BILDER



UNI

de	15 mm-Konnektor (UNI)
en	15 mm-connector (UNI)
fr	Connecteur 15 mm (UNI)
nl	15 mm-connector (UNI)
es	Conector de 15 mm (UNI)
sv	15 mm anslutning (UNI)



PHON

de	Sprechventil (PHON)
en	Speaking valve (PHON)
fr	Valve de phonation (PHON)
nl	Spreekventiel (PHON)
es	Válvula fonatoria (PHON)
sv	Talventil (PHON)



VARIO

de	15 mm-Drehkonnektor (VARIO)
en	15 mm swivel connector (VARIO)
fr	Connecteur rotatif 15 mm (VARIO)
nl	15 mm swivel connector (VARIO)
es	Conector rotatorio de 15 mm (VARIO)
sv	15 mm vridkoppling (VARIO) (VARIO)



CUFF

de	Mit CUFF
en	With CUFF
fr	Avec ballonnet (CUFF)
nl	Met CUFF
es	Con balón (CUFF)
sv	Med CUFF



KOMBI

de	22 mm-Konnektor (KOMBI)
en	22 mm-connector (KOMBI)
fr	Connecteur 22 mm (KOMBI)
nl	22 mm-connector (KOMBI)
es	Conector de 22 mm (KOMBI)
sv	22 mm anslutning (KOMBI)



SUCTION

de	Mit Absaugvorrichtung (SUCTION)
en	With suction line (SUCTION)
fr	Avec dispositif d'aspiration (SUCTION)
nl	Met afzuigkanaal (SUCTION) voor subglottisch afzuigen van secreet boven de cuff
es	Con línea de succión (SUCTION)
sv	Med sugslang (SUCTION)



LINGO

de	Fensterung oder Siebung (LINGO)
en	Fenestration or sieving (LINGO)
fr	Fenêtrage ou perforations (LINGO)
nl	Vensteren of zeven (LINGO)
es	Fenestration o perforación (LINGO)
sv	Fönsteröppning eller siktning (LINGO)

O₂

de	Mit einem O ₂ -Anschluss
en	With oxygen support
fr	Avec raccord O ₂
nl	Met zuurstofaansluiting (O ₂ -support)
es	Con puerto de oxígeno
sv	Med O ₂ -anslutning



FLAT-PHON

de	Flachprofil mit Sprechventil (FLAT-PHON)
en	Flat profile with speaking valve (FLAT-PHON)
fr	Profilé plat avec valve de phonation (FLAT-PHON)
nl	Vlak profiel met spreekventiel (FLAT-PHON)
es	Perfil plano con válvula fonatoria (FLAT-PHON)
sv	Platt profil med talventil (FLAT-PHON)



MRT

de	MRT geeignet
en	MRT suitable
fr	IRM-compatible
nl	Geschikt voor MRT
es	Compatible con resonancia Magnética
sv	MR-kompatibel

de Falls zutreffend, sind diese Piktogramme auf dem Datenblatt zu finden.
 en If applicable, these pictograms can be found on the data sheet.
 fr Le cas échéant, ces pictogrammes se trouvent sur la fiche technique.

nl Indien van toepassing zijn deze pictogrammen te vinden op het gegevensblad.
 es Si procede, estos pictogramas se encuentran en la ficha técnica.
 sv I förekommande fall finns dessa piktogram på databladet.