

MUCOPROTECT® CUFF PRESSURE GAUGE



- ▶ For cuff pressure monitoring of tracheostomy tubes and endotracheal tubes with cuff
- ▶ Large surface allows precise adjustment to ensure maximum control accuracy
- ▶ Optimal pressure range (20-30 cm H₂O) highlighted in blue colour
- ▶ Lateral pressure release button
- ▶ Can be used for all control balloons with Luer-Lock-connector (market standard)
- ▶ Connection tube (Length: 75 cm) included
- ▶ Connection tube (**MUCOPROTECT® Connection tube** / REF 19506) also separately available

SCOPE OF DELIVERY

**MUCOPROTECT®
CUFF PRESSURE GAUGE** REF 19500

1 x Cuff pressure gauge
1 x Connection tube
1 x Instruction for use

Optional available:
**MUCOPROTECT®
CONNECTION TUBE (75 cm)** REF 19506

CE 0482