

LARYVOX® ADAPT

REF 48920



- ▶ Permits the connection of aids such as HMEs with a 22-mm connector (**KOMBI**) to a tracheostomy tube with a 15-mm connector (**UNI**)
- ▶ There is at one end a 22-mm connector, which is used for attaching an HME and has an adapter at the other end which can take a plug connection with a 15-mm connector (**UNI**)

SCOPE OF DELIVERY

LARYVOX® ADAPT

REF 48920

3 x adaptor, individually packed

1 x Instruction for use

