

VOICETEC®

REF 77350



- ▶ Enables excellent amplification of a low esophageal voice
- ▶ For a low voice due to impaired function of the vocal cords
- ▶ Ideal for amplifying your voice in noisy environments
- ▶ For use as an external speaker for sound sources such as smartphones, tablets, computers, etc.
- ▶ Clear sound quality
- ▶ Continuous volume control
- ▶ Long operating time up to 22 hours
- ▶ Wearable with adjustable strap
- ▶ Low weight (125 g)
- ▶ Dimensions: 9 x 7 x 2,5 cm
- ▶ Bluetooth enabled

SCOPE OF DELIVERY**VOICETEC®****REF 77350**

1x Voice Amplifier
 1x Cable-bound headset microphone
 1x Waist band
 1x USB Typ C-charging cable
 1x AUX audio cable
 1x Instructions for use

optionally available

**Headset Microphone
VOICETEC®****REF 77350-03**