

ABSORBELLIN®

sterile

REF 30850



REF 30851



REF 30852



REF 30853



REF 30854



- ▶ High-quality tracheal compresses, with anatomically adapted rounding for optimum comfort
- ▶ All **Absorbellin®** compresses are made of hydrophilic polyurethane foam, which reliably absorbs liquids
- ▶ **Absorbellin®** compresses lie softly on sensitive skin and cushion the sensitive skin areas thanks to the dense polyurethane foam
- ▶ **Absorbellin® PED** is specially made for newborns and infants
- ▶ **Absorbellin® Drain** for cushioning drains, probes, catheters and for absorbing wound fluids
- ▶ All compresses are sterile
- ▶ **SLIT** variants with slit for changing tracheostomy tubes or drains in place

SCOPE OF DELIVERY	
Absorbellin®, slit, sterile 10 x Tracheal compresses	REF 30850
or Absorbellin®, sterile 10 x Tracheal compresses	REF 30851
or Absorbellin® PED, sterile 10 x Tracheal compresses	REF 30852
or Absorbellin® PED, slit, sterile 10 x Tracheal compresses	REF 30853
or Absorbellin® Drain, slit, sterile 10 x Drainage compresses	REF 30854

Product	Size	PU	REF
Absorbellin®, slit, sterile	9 x 10 cm	10	30850
Absorbellin®, sterile	9 x 10 cm	10	30851
Absorbellin® PED, sterile	6 x 6,5 cm	10	30852
Absorbellin® PED, slit, sterile	6 x 6,5 cm	10	30853
Absorbellin® Drain, slit, sterile	7,5 x 9 cm	10	30854

CE 0482