

NECKFIX® ENDO

REF 32614



REF 32616



REF 32613



REF 32615



- ▶ Fixation strap for endotracheal tubes, enables secure fixation of endotracheal tubes
- ▶ Lightweight foam core, encased in a skin-friendly velour layer
- ▶ The tube holder is approx. 6 mm thick and 15 mm wide and is available in different versions:
 - **NECKFIX® ENDO 72 cm**
 - **NECKFIX® ENDO SAFETY K 72 cm** – with additional safety adhesive Velcro strip
 - **NECKFIX® ENDO 85 cm**
 - **NECKFIX® ENDO SAFETY K 85 cm** – with additional safety adhesive Velcro strip
- ▶ Enables quick adjustment to the desired neck size by means of a Velcro fastener
- ▶ Available in various lengths
- ▶ Colour: white-turquoise

SCOPE OF DELIVERY

NECKFIX® ENDO 72 cm REF 32614
10 x 10 tube holder

or

NECKFIX® ENDO Safety K 72 cm REF 32616
10 x 10 tube holder with additional safety adhesive Velcro strip

or

NECKFIX® ENDO 85 cm REF 32613
10 x 10 tube holder

or

NECKFIX® ENDO Safety K 85 cm REF 32615
10 x 10 tube holder with additional safety adhesive Velcro strip

Optionally available:

NECKFIX® Endo Velcro Strips REF 32617
100 x adhesive Velcro strips

