

FAHL® HAND FIX

REF 32470



REF 32472



REF 32482



REF 32485



- ▶ Padded cuff with Velcro fastener and plastic ring for hand or foot fixation
- ▶ Soft cuff with skin-friendly velour coating for comfortable wear and a secure fit
- ▶ For self-securing patients
- ▶ The cuffs are approx. 6 mm thick and available in different versions:
 - **NEO (REF 32470):**
3 x 12 cm (W x L), strap length approx. 2 x 65 cm
 - **PED (REF 32472):**
4 x 17 cm (W x L), strap length approx. 2 x 65 cm
 - **Normal (REF 32482):**
7 x 25 cm (W x L), strap length approx. 2 x 85 cm
 - **Large (L) (REF 32485):**
7 x 30 cm (W x L), strap length approx. 2 x 85 cm
- ▶ Colour: white-turquoise / white

SCOPE OF DELIVERY

FAHL® Hand Fix NEO, 3 x 12 cm

15 x 2 Fixation cuffs

REF 32470

or

FAHL® Hand Fix PED, 4 x 17 cm

15 x 2 Fixation cuffs

REF 32472

or

FAHL® Hand Fix, 7 x 25 cm

8 x 2 Fixation cuffs

REF 32482

or

FAHL® Hand Fix L, 7 x 30 cm

8 x 2 Fixation cuffs

REF 32485

