


TRACHEOTEC® VARIO CUFF SUCTION

REF 18930



- ▶ Tracheostomy tube with low-pressure cuff (**CUFF**) and integrated suction tube (**SUCTION**) to remove secretions above the **CUFF**
- ▶ The suction tube is integrated into the inner wall of the cannula tube to avoid additional irritation of the trachea or tracheostoma
- ▶ Soft, flexible cannula tag with imprint of the essential specifications of the cannula
- ▶ With 15 mm swivel connector (**VARIO**) to buffer tensile or torsional forces that may be generated by the respiratory tube
- ▶ With X-ray contrast strip
- ▶ Set includes tube holder, an obturator and cannula passport

SCOPE OF DELIVERY

TRACHEOTEC® VARIO CUFF SUCTION

REF 18930-xx

- 1 x Tracheostomy tube with **CUFF** and suction line (**SUCTION**)
- 1 x Tube holder
- 1 x Obturator
- 1 x Cannula ID card
- 1 x Instructions for Use

Please state the size after the article number when ordering! You will find the size chart in the appendix.



VARIO



CUFF




SUCTION

CE 0482

TRACHEOTEC® VARIO CUFF SUCTION

TECHNICAL SKETCH



TRACHEOTEC® VARIO CUFF SUCTION

REF 18930-xx

Size	OD1	ID1	CD	NL	OL
03	4,7	3,0	15,0	43,3	47,2
035	5,3	3,5	15,0	46,2	50,6
04	6,0	4,0	19,0	52,1	57,1
045	6,7	4,5	19,0	52,1	57,7
05	7,3	5,0	19,0	55,1	61,1
055	8,0	5,5	19,0	57,0	63,6
06	8,7	6,0	22,0	59,9	67,1
065	9,3	6,5	22,0	65,8	73,5
07	10,0	7,0	30,0	68,7	77,0
075	10,7	7,5	30,0	72,6	81,5
08	11,3	8,0	30,0	78,5	87,9
085	12,0	8,5	30,0	80,5	90,4
09	12,7	9,0	30,0	84,4	94,9
10	14,0	10,0	30,0	93,7	105,3

Legend

All specifications in mm

OD1 = Outer diameter at the neck flange

ID1 = Inner diameter at the cannula tip of the outer cannula

CD = CUFF diameter

NL = Length over the central arch

OL = Length over the outer curve

PIKTOGRAMME / PICTOGRAMS / PICTOGRAMMES / PICTOGRAMMEN / PICTOGRAMAS / BILDER



UNI

de	15 mm-Konnektor (UNI)
en	15 mm-connector (UNI)
fr	Connecteur 15 mm (UNI)
nl	15 mm-connector (UNI)
es	Conector de 15 mm (UNI)
sv	15 mm anslutning (UNI)



PHON

de	Sprechventil (PHON)
en	Speaking valve (PHON)
fr	Valve de phonation (PHON)
nl	Spreekventiel (PHON)
es	Válvula fonatoria (PHON)
sv	Talventil (PHON)



VARIO

de	15 mm-Drehkonnektor (VARIO)
en	15 mm swivel connector (VARIO)
fr	Connecteur rotatif 15 mm (VARIO)
nl	15 mm swivel connector (VARIO)
es	Conector rotatorio de 15 mm (VARIO)
sv	Ytterkanyl inklusive 15 mm vridanslutning (VARIO)



CUFF

de	Mit CUFF
en	With CUFF
fr	Avec ballonnet (CUFF)
nl	Met CUFF
es	Con balón (CUFF)
sv	Med CUFF



KOMBI

de	22 mm-Konnektor (KOMBI)
en	22 mm-connector (KOMBI)
fr	Connecteur 22 mm (KOMBI)
nl	22 mm-connector (KOMBI)
es	Conector de 22 mm (KOMBI)
sv	22 mm anslutning (KOMBI)



SUCTION

de	Mit Absaugvorrichtung (SUCTION)
en	With suction line (SUCTION)
fr	Avec dispositif d'aspiration (SUCTION)
nl	Met afzuigkanaal (SUCTION) voor subglottisch afzuigen van secreet boven de cuff
es	Con línea de succión (SUCTION)
sv	Med sugslang (SUCTION)



LINGO

de	Fensterung oder Siebung (LINGO)
en	Fenestration or sieving (LINGO)
fr	Fenêtrage ou perforations (LINGO)
nl	Vensteren of zeven (LINGO)
es	Fenestration o perforación (LINGO)
sv	Fönsteröppning eller siktning (LINGO)

O₂

de	Mit einem O ₂ -Anschluss
en	With oxygen support
fr	Avec raccord O ₂
nl	Met zuurstofaansluiting (O ₂ -support)
es	Con puerto de oxígeno
sv	Med O ₂ -anslutning



FLAT-PHON

de	Flachprofil mit Sprechventil (FLAT-PHON)
en	Flat profile with speaking valve (FLAT-PHON)
fr	Profilé plat avec valve de phonation (FLAT-PHON)
nl	Vlak profiel met spreekventiel (FLAT-PHON)
es	Perfil plano con válvula fonatoria (FLAT-PHON)
sv	Platt profil med talventil (FLAT-PHON)



MRT

de	MRT geeignet
en	MRT suitable
fr	IRM-compatible
nl	Geschikt voor MRT
es	Compatible con resonancia Magnética
sv	MR-kompatibel

de Falls zutreffend, sind diese Piktogramme auf dem Datenblatt zu finden.
 en If applicable, these pictograms can be found on the data sheet.
 fr Le cas échéant, ces pictogrammes se trouvent sur la fiche technique.

nl Indien van toepassing zijn deze pictogrammen te vinden op het gegevensblad.
 es Si procede, estos pictogramas se encuentran en la ficha técnica.
 sv I förekommande fall finns dessa piktogram på databladet.