

COMBIPHON® BLUE/ COMBIPHON® O₂ BLUE

REF 27131-02



CE

REF 27132-02



CE 0482

- ▶ Speaking valve - low weight
- ▶ The blue colour distinguishes it unmistakably from other attachments, such as a HME
- ▶ With silicone membrane
- ▶ The valve membrane is closed in the basic position and opens briefly during inhalation. Advantage: a high blowing pressure is not required to close the valve and thus enable speaking.
- ▶ Reduced fatigue for the patient and improved speaking quality
- ▶ With multi-adapter: fits a 15 mm connector (**UNI**) and also a 22 mm adapter (**KOMBI**)
- ▶ For attachment to a tracheostomy tube or base plate
- ▶ Only for tracheostomees, not for laryngectomees
- ▶ The speaking valve is intended for use only and exclusively in combination with a phonation cannula.

SCOPE OF DELIVERY

COMBIPHON® BLUE**REF 27131-02**

1 x speaking valve
1 x Instructions for Use

or

COMBIPHON® O₂ BLUE**REF 27132-02**

1 x speaking valve
1 x Instructions for Use



UNI



KOMBI

O₂

The COMBIPHON® speaking valve must only be used by patients with partially or fully preserved larynx!

Product	PU	REF
COMBIPHON® BLUE	1 piece	27131-02
COMBIPHON® O ₂ BLUE	1 piece	27132-02