

LARYVOX® HME

REF 49800



CE

REF 49801




CE 0482

REF 49802



CE 0482



- ▶ Open pore polyurethane foam HME filter
- ▶ Beige colour to provide inconspicuous design
- ▶ Lightweight; low profile
- ▶ Stable housing – can be closed by light finger occlusion to enable laryngectomized patients to speak
- ▶ With multi-adapter to fit on 15 mm-connectors (**UNI**) and on 22 mm adapters (**KOMBI**) to connect to any conventional base plates (e.g. **LARYVOX® TAPE**) or tracheostomy tubes with 15 mm-connectors or silicone tubes for laryngectomees (e.g. **LARYNGOTEC®**)
- ▶ Available in three variants:
 - **REF 49800**: HME filter
 - **REF 49801**: O₂-Connection
 - **REF 49802**: HME filter with O₂-Connection

SCOPE OF DELIVERY

LARYVOX® HME **REF 49800**

6 x 5 Filter cassettes
1 x Instruction for use

or

LARYVOX® O₂ Connect **REF 49801**

5 x O₂-connector, individually packed
1 x Instruction for use

or

LARYVOX® O₂ HME **REF 49802**

6 x 5 Filter cassettes
1 x Instruction for use