

# LARYNGOTEC® KOMBI STOMA BUTTON FIX

REF 115950



- ▶ Stoma button to keep the tracheostoma of laryngectomized patients open while minimizing the contact surface with the delicate tracheal mucosa
- ▶ For long-term use, can be used up to 3-6 months
- ▶ Medical grade silicone: soft, pliable, lightweight
- ▶ With 22 mm adapter (**KOMBI**) to connect to HMEs and/or speaking valves
- ▶ With two lateral eyelets to fasten a tube holder to prevent dislocation of the stoma button
- ▶ In two different lengths
- ▶ Transparent, inconspicuous design

## SCOPE OF DELIVERY

### LARYNGOTEC® KOMBI STOMA BUTTON FIX

Size 8-12

REF 115950

- 1 x Stoma button
- 1 x Including tube holder (OPTIFLAUSCH®)
- 1 x product ID card
- 1 x Instruction for use

or

### LARYNGOTEC® KOMBI STOMA BUTTON

FIX SHORT, Size 8-12

REF 115960

- 1 x Stoma button
- 1 x Including tube holder (OPTIFLAUSCH®)
- 1 x product ID card
- 1 x Instruction for use

When ordering, please state the size after the item number!



KOMBI




MRT

CE 0482

## LARYNGOTEC® KOMBI STOMA BUTTON FIX

## Technical Data



## LARYNGOTEC® KOMBI STOMA BUTTON FIX REF 115950-xx

| Size | OD   | ID   | NL   | REF       |
|------|------|------|------|-----------|
| 8    | 11,0 | 8,0  | 22,0 | 115950-08 |
| 9    | 12,0 | 9,0  | 22,0 | 115950-09 |
| 10   | 13,0 | 10,0 | 22,0 | 115950-10 |
| 11   | 14,0 | 11,0 | 22,0 | 115950-11 |
| 12   | 15,0 | 12,0 | 22,0 | 115950-12 |

## LARYNGOTEC® KOMBI STOMA BUTTON FIX SHORT REF 115960-xx

| Size | OD   | ID   | NL   | REF       |
|------|------|------|------|-----------|
| 8    | 11,0 | 8,0  | 15,0 | 115960-08 |
| 9    | 12,0 | 9,0  | 15,0 | 115960-09 |
| 10   | 13,0 | 10,0 | 15,0 | 115960-10 |
| 11   | 14,0 | 11,0 | 15,0 | 115960-11 |
| 12   | 15,0 | 12,0 | 15,0 | 115960-12 |

## Legend

All lengths in mm

OD = Outer diameter at the neck flange

ID = Inner diameter at the canula tip of the outer cannula

NL = Length over the middle curve

# PIKTOGRAMME / PICTOGRAMS / PICTOGRAMMES / PICTOGRAMMEN / PICTOGRAMAS / BILDER



UNI

|    |                         |
|----|-------------------------|
| de | 15 mm-Konnektor (UNI)   |
| en | 15 mm-connector (UNI)   |
| fr | Connecteur 15 mm (UNI)  |
| nl | 15 mm-connector (UNI)   |
| es | Conector de 15 mm (UNI) |
| sv | 15 mm anslutning (UNI)  |



VARIO

|    |                                                   |
|----|---------------------------------------------------|
| de | 15 mm-Drehkonnektor (VARIO)                       |
| en | 15 mm-swivel connector (VARIO)                    |
| fr | Connecteur rotatif 15 mm (VARIO)                  |
| nl | 15 mm-swivel-connector (VARIO)                    |
| es | Conector rotatorio de 15 mm (VARIO)               |
| sv | Ytterkanyl inklusive 15 mm vridanslutning (VARIO) |



KOMBI

|    |                           |
|----|---------------------------|
| de | 22 mm-Konnektor (KOMBI)   |
| en | 22 mm-connector (KOMBI)   |
| fr | Connecteur 22 mm (KOMBI)  |
| nl | 22 mm-connector (KOMBI)   |
| es | Conector de 22 mm (KOMBI) |
| sv | 22 mm anslutning (KOMBI)  |



LINGO

|    |                                       |
|----|---------------------------------------|
| de | Fensterung oder Siebung (LINGO)       |
| en | Fenestration or sieving (LINGO)       |
| fr | Fenêtrage ou perforations (LINGO)     |
| nl | Vensteren of zeven (LINGO)            |
| es | Fenestration o perforación (LINGO)    |
| sv | Fönsteröppning eller siktning (LINGO) |




FLAT-PHON

|    |                                                  |
|----|--------------------------------------------------|
| de | Flachprofil mit Sprechventil (FLAT-PHON)         |
| en | Flat profile with speaking valve (FLAT-PHON)     |
| fr | Profilé plat avec valve de phonation (FLAT-PHON) |
| nl | Vlak profiel met spreekventiel (FLAT-PHON)       |
| es | Perfil plano con válvula fonatoria (FLAT-PHON)   |
| sv | Platt profil med talventil (FLAT-PHON)           |



PHON

|    |                           |
|----|---------------------------|
| de | Sprechventil (PHON)       |
| en | Speaking valve (PHON)     |
| fr | Valve de phonation (PHON) |
| nl | Spreekventiel (PHON)      |
| es | Válvula fonatoria (PHON)  |
| sv | Talventil (PHON)          |



CUFF

|    |                       |
|----|-----------------------|
| de | Mit CUFF              |
| en | With CUFF             |
| fr | Avec ballonnet (CUFF) |
| nl | Met CUFF              |
| es | Con balón (CUFF)      |
| sv | Med CUFF              |



SUCTION

|    |                                                                                 |
|----|---------------------------------------------------------------------------------|
| de | Mit Absaugvorrichtung (SUCTION)                                                 |
| en | With suction line (SUCTION)                                                     |
| fr | Avec dispositif d'aspiration (SUCTION)                                          |
| nl | Met afzuigkanaal (SUCTION) voor subglottisch afzuigen van secreet boven de cuff |
| es | Con línea de succión (SUCTION)                                                  |
| sv | Med sugslang (SUCTION)                                                          |

O<sub>2</sub>

|    |                                                   |
|----|---------------------------------------------------|
| de | Mit einem O <sub>2</sub> -Anschluss               |
| en | With oxygen support                               |
| fr | Avec raccord O <sub>2</sub>                       |
| nl | Met zuurstofaansluiting (O <sub>2</sub> -support) |
| es | Con puerto de oxígeno                             |
| sv | Med O <sub>2</sub> -anslutning                    |



MRT

|    |                                     |
|----|-------------------------------------|
| de | MRT geeignet                        |
| en | MRT suitable                        |
| fr | IRM-compatible                      |
| nl | Geschikt voor MRT                   |
| es | Compatible con resonancia Magnética |
| sv | MR-kompatibel                       |

de Falls zutreffend, sind diese Piktogramme auf dem Datenblatt zu finden.  
 en If applicable, these pictograms can be found on the data sheet.  
 fr Le cas échéant, ces pictogrammes se trouvent sur la fiche technique.

nl Indien van toepassing zijn deze pictogrammen te vinden op het gegevensblad.  
 es Si procede, estos pictogramas se encuentran en la ficha técnica.  
 sv I förekommande fall finns dessa piktogram på databladet.