

DYSSUCTION® SUCTION DEVICE M

REF 67900



REF 100705-2 (optional)



- ▶ Portable suction device for subglottic suction
- ▶ Disposable secretion container with integrated suction tube for hygienic and germ-reduced container changes
- ▶ The disposable secretion container can be used for up to 7 days
- ▶ Integrated activated charcoal filter for odour reduction
- ▶ Integrated gelling agent for binding secretions
- ▶ Digital display in the languages German, English, French, Spanish, Italian, Dutch
- ▶ Individually adjustable
 - Vacuum: -60 mbar to -300 mbar (in 10 mbar increments)
 - Suction time: 10 - 60 sec (setting in 1-second increments)
 - Interval duration: 3 - 60 min (setting in 1-minute increments)

SCOPE OF DELIVERY**DYSSUCTION® SUCTION DEVICE M REF 67900**

- 1 x suction device incl. disposable secretion container (250 ml) with suction tube with stepped connector
- 1 x disposable secretion container (250 ml) with suction tube with stepped connector
- 1 x power supply unit **REF 67912**
- 1 x Instructions for Use - Quick Guide
- 1 x Instructions for Use

Optionally available:

- 30 x disposable secretion containers (250 ml) with stepped connector **REF 67910**
- 1 x holder variable DYSSuction Suction Device M **REF 100705-2**

FAHL

DYSSUCTION® SUCTION DEVICE M



DEVICE OVERVIEW

- ① Control panel
- ② Removable secretion container
- ③ Suction tube
- ④ Charging connection

TECHNICAL DATA

Vacuum:	-60 mbar to -300 mbar (in 10 mbar increments)
Suction time:	10 - 60 sec (in 1-second increments)
Interval duration:	3 - 60 min (in 1-minute increments)
Unit suction capacity	max. 8 l/min
Dimensions basic device (H x W x D)	165 x 220 x 90 mm
Weight	1.2 kg
Operating time	Mains operation: permanent operation Battery operation: approx. 18 h depending on the load on the motor
Charging time	6 – 7 h with empty battery
Sound emission	35 dB (A)